**Monoclonal Mouse  
Anti-Human  
HLA-ABC Antigen/RPE  
Clone W6/32  
Code No./ Code/ Code-Nr . R7000**

ENGLISH							
<b>Intended use</b>	For in vitro diagnostic use. R7000 is intended for use in flow cytometry. Anti-HLA-ABC antigen is relevant for the study of HLA class I expression in patients with solid tumors (1-4). R7000 is not intended for use in tissue typing. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.						
<b>Introduction</b>	The human leucocyte antigen (HLA) system, originally discovered as the result of a transfusion reaction, is now known to play a crucial role in many areas of clinical medicine. The HLA molecules are encoded by a cluster of tightly linked genes located on the short arm of chromosome 6. Based on some of the structural and functional characteristics of the genes, the region has been divided into three: HLA class I, Class II and class III regions. The class I region contains genes encoding for the heavy chain of the HLA class I molecules. The HLA class I genes have been classified according to their structure, expression and function as classical (HLA-A, B and C) and non-classical (HLA-E, F and G). Both classical and non-classical HLA class I genes encode a heavy $\alpha$ -chain, of approximately 43 kDa, non-covalently linked to a non-polymorphic light chain, the $\beta$ 2-microglobulin which is encoded by a gene on chromosome 15 (5). The main function of the HLA-A, B and C molecules is to present antigenic peptides, derived primarily but not exclusively from endogenous proteins, to CD8+ T-cells. HLA molecules are also known to be associated with a variety of autoimmune, non-autoimmune and infectious diseases and to restrict the antibody response to certain antigens and vaccines (5). HLA-A, -B and -C antigens are widely distributed on most human nucleated cells. However, the intensity of expression varies considerably, some cells being only weakly positive, e.g. thyroid and muscle cells, and others negative, e.g. cells of the exocrine pancreas and villous trophoblast cells (6). The intensity of HLA-ABC antigens may also be altered in pathological states. It has been described that malignant cells may lose HLA-ABC (1-4), whereas hepatocytes in alcoholic hepatitis, biliary cirrhosis and chronic active hepatitis express HLA-ABC, in contrast to normal liver hepatocytes on which HLA-ABC was not detected (7).						
<b>Reagent provided</b>	The Anti-HLA-ABC Antigen conjugate, R7000, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 100 tests (10 $\mu$ L of conjugate for up to 10 <sup>6</sup> leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG2a, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.						
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>R7000</td> <td>RPE (R-Phycoerythrin)</td> <td>X0950</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	R7000	RPE (R-Phycoerythrin)	X0950
Antibody Code No.	Fluorochrome	Negative Control Code No.					
R7000	RPE (R-Phycoerythrin)	X0950					
<b>Immunogen</b>	Cell membranes of human tonsil cells (8).						
<b>Specificity</b>	Anti-HLA-ABC Antigen, W6/32, is directed against a monomorphic epitope on the 45 kDa polypeptide products of the HLA-A, B and C loci (8). Anti-HLA-ABC Antigen, W6/32, reacts with all nucleated cells in peripheral blood or tonsil cell preparations. These include polymorphs, monocytes, lymphocytes and eosinophils. The antibody does not react with erythrocytes (8). Flow cytometric analysis of samples from epithelial ovarian cancer patients demonstrated that Anti-HLA-ABC Antigen, W6/32, labelling intensity may be used as a prognostic indicator. Low HLA-ABC antigen expression may be a poor prognostic indicator among patients with aneuploid ovarian cancers (3).						
<b>Precautions</b>	<ol style="list-style-type: none"> <li>For professional users.</li> <li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>As with any product derived from biological sources, proper handling procedures should be used.</li> </ol>						
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.						
<b>Staining procedure</b>	<ol style="list-style-type: none"> <li>Collect venous blood into a test tube containing an anticoagulant.</li> <li>Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.</li> <li>Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.</li> <li>Mix 100 <math>\mu</math>L cell suspension with 10 <math>\mu</math>L R7000.</li> <li>Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).</li> <li>Incubate in the dark at 4 °C for 30 minutes.</li> <li>Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.</li> <li>Analyse on a flow cytometer.</li> </ol> <p>It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p>						
(104391-002)	R7000/EFG/MBH/04.03.05 p. 1/4						

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

1. Bander NH, Yao D, Liu H, Chen YT, Steiner M, Zuccaro W, Moy P. MHC class I and II expression in prostate carcinoma and modulation by interferon-alpha and -gamma. Prostate 1997;33:233-9.
2. Rockett JC, Damton SJ, Crocker J, Matthews HR, Morris AG. Expression of HLA-ABC, HLA-DR and intercellular adhesion molecule-1 in oesophageal carcinoma. J Clin Pathol 1995;48:539-44.
3. Moore DH, Fowler WC, Olafsson K. Class I histocompatibility antigen expression: a prognostic factor for aneuploid ovarian cancers. Gynecol Oncol 1990;38:458-61.
4. Fleming KA, McMichael A, Morton JA, Woods J, McGee JO. Distribution of HLA class I antigens in normal human tissue and mammary cancer. J Clin Pathol 1981;34:779-84.
5. Navarrete CV. The HLA system in blood transfusion. Baillière's Clin Haematol 2000;13:511-32.
6. Daar AS, Fuggle SV, Fabre JW, Ting A, Morris PJ. The detailed distribution of HLA-A, B, C antigens in normal human organs. Transplantation 1984;38:287-92.
7. Barbatis C, Woods J, Morton JA, Fleming KA, McMichael A, McGee JO. Immunohistochemical analysis of HLA (A,B,C) antigens in liver disease using a monoclonal antibody. Gut 1981;22:985-91.
8. Barnstable CJ, Bodmer WF, Brown G, Galfre G, Milstein C, Williams AF, et al. Production of monoclonal antibodies to group A erythrocytes, HLA and other human cell surface antigens - new tools for genetic analysis. Cell 1978;14:9-20.

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