**Monoclonal Mouse Anti-Human CD22/APC, Clone 4KB128**  
**Monoclonal Mouse Anti-Human CD22/FITC, Clone 4KB128**  
**Monoclonal Mouse Anti-Human CD22/RPE, Clone 4KB128**

**C7281**  
**F7060**  
**R7061**

## ENGLISH

**Intended use** For in vitro diagnostic use.  
 C7281, F7060 and R7061 are intended for use in flow cytometry. CD22 is a pan-B cell reagent enabling detection of normal and neoplastic B cells in peripheral blood. CD22 is a recommended marker for B-cell lineage acute leukaemia. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Synonyms for antigen** B lymphocyte cell adhesion molecule (BL-CAM), Lyb-8.2 (1, 2).

**Introduction** CD22 is a 140 kDa single chain type 1 transmembrane glycoprotein (3). Cloning of the CD22 gene has revealed that CD22 is a member of the IgG superfamily and homologous to several other proteins, including myelin basic protein and members of the carcinoembryonic antigen (CEA) family (2). The structural features of CD22 suggest that CD22 would mediate intercellular adhesion, however, further work is needed to understand the in vivo function of CD22 and to identify physiologically significant ligands (1).

CD22 expression is restricted to normal and neoplastic B cells and is absent from other haemopoietic cell types (4). CD22 is first expressed in pro-B and pre-B cells, largely as a cytoplasmic protein, and only later at the mature B-cell stage is it found on the cell surface (1,4). The antigen is lost during the terminal stages of differentiation prior to the plasma cell stage. Most precursor-B-cell acute lymphoblastic leukaemias (ALL) express both cytoplasmic and membrane forms of CD22, thus prominent surface expression of CD22 was observed in 90% of 52 precursor-B-ALL (4).

**Reagent provided** The Anti-CD22 conjugates, C7281, F7060 and R7061, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10<sup>6</sup> leucocytes from normal human peripheral blood).

**Isotype:** IgG1, kappa. **Conjugate concentration mg/L:** See label on vial.

| Antibody Code No. | Fluorochrome                               | Control Reagent Code No. |
|-------------------|--|--------------------------|
| C7281             | APC (Allophycocyanin)                      | X0968                    |
| F7060             | FITC (Fluorescein Isothiocyanate Isomer 1) | X0927                    |
| R7061             | RPE (R-Phycoerythrin)                      | X0928                    |

**Immunogen** Neoplastic cells from hairy cell leukaemia (5).

**Specificity** Anti-CD22, 4KB128, was included in the Third and Fifth International Workshops and Conferences on Human Leucocyte Differentiation Antigens (4, 5), and studies by different laboratories confirmed its reactivity with the CD22 antigen.

Anti-CD22, 4KB128, was shown to detect surface membrane CD22 on precursor B cells (4). A weak staining of neutrophils in blood smears may occasionally be observed with the antibody.

Similarly, a weak staining of neutrophils may be seen when F7060 is used in flow cytometry for intracellular staining of leucocytes from normal human peripheral blood.

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage** Store in the dark at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

1. Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
2. Add 10 µL of fluorochrome-conjugated anti-CD22 and mix gently by using a vortex mixer.
3. Incubate the tube in the dark at 2–8 °C for 30 minutes or at room temperature (20–25 °C) for 15–30 minutes.
4. Add 100 µL of Uti-Lyse™ Reagent A (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 1 mL of Uti-Lyse™ Reagent B (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
7. Add 2 mL of PBS (Dako code No. S3024) to the tube and resuspend the cells by using a vortex mixer.
8. Repeat step 6.
9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.
10. Analyse the sample on a flow cytometer or store it in the dark at 2–8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

**Procedural notes**

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the isotype and fluorochrome of the conjugated antibody. Recommended control reagents are shown in the table above.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as Dako EasyLyse™, code No. S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis.

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

1. Tedder TF, Tuscano J, Sato S, Kehrl JH. CD22, a B lymphocyte-specific adhesion molecule that regulates antigen receptor signaling. *Annu Rev Immunol* 1997;15:481-504.
2. Wilson GL, Fox CH, Fauci AS, Kehrl JH. cDNA cloning of the B cell membrane protein CD22: a mediator of B-B cell Interactions. *J Exp Med* 1991;173:137-46.
3. Engel P. CD guide. CD22. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. *Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 767-8.*
4. Kehrl JH. B6. CD22 workshop panel report. In: Schlossman SF, Bounsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press. Volume 1. 1995. p. 523-5.*
5. Moldenhauer G, Schwartz R, Dörken B, Hämmerling GJ. B2.2. Biochemical characterization and epitope analysis of B-lymphocyte-specific surface antigens defined by clustering Workshop monoclonal antibodies. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. *Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 378-82.*

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Agilent Technologies Singapore (International) Pte Ltd.