**Monoclonal Mouse  
Anti-Human  
CD45RA/RPE  
Clone 4KB5  
Code No./ Code/ Code-Nr. R 7086**

**ENGLISH**

<b>Intended use</b>	For in vitro diagnostic use. R 7086 is intended for use in flow cytometry. Anti-CD45RA, 4KB5, reacts with most B cells in peripheral blood. It also recognizes a small subpopulation of T cells and monocytes (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
<b>Introduction</b>	CD45 is a transmembrane glycoprotein expressed on most nucleated cells of haematopoietic origin. CD45, encoded by a single gene mapped to chromosome 1, has various isoforms, based on differential splicing of exons 4, 5 and 6. On human leucocytes, five different isoforms of CD45, named ABC, AB, BC, B and 0, have been identified. These isoforms are recognized by CD45RA, CD45RB, CD45RC and CD45R0 antibodies. Antibodies to CD45RA recognize the extracellular A domain in the two isoforms, ABC (Mr 220 000), and AB (Mr 200 000). All the CD45 isoforms share the same intracellular segment, which has been shown to have tyrosine phosphatase activity. Various leucocytes express characteristic CD45 isoforms, and T cells differentially express CD45 isoforms at various stages of their development and activation. B cells predominantly express the ABC isoform, and monocytes and dendritic cells predominantly express the B and 0 isoforms. Granulocytes principally express only the B and 0 isoforms (2). Discrimination between naive and memory T cells has been made by establishing their CD45 isoform. It is believed that CD45R0 is expressed preferentially on memory T cells and CD45RA on naive T cells (3).
<b>Reagent provided</b>	The Anti-CD45RA conjugate, R 7086, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 <sup>6</sup> leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.
<b>Immunogen</b>	Human hairy cell leukaemia cells (4).
<b>Specificity</b>	The antibody was included at the Third (4) and clustered as anti-CD45RA at the Fourth (5) and Fifth (6) International Workshops and Conferences on Human Leucocyte Differentiation Antigens.
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
<b>Staining procedure</b>	<ol style="list-style-type: none"> <li>1. Collect venous blood into a test tube containing an anticoagulant.</li> <li>2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.</li> <li>3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.</li> <li>4. Mix 100 µL cell suspension with 10 µL RPE-conjugated Anti-CD45RA.</li> <li>5. Mix 100 µL cell suspension with 10 µL of the appropriate negative control. The recommended negative control reagent for R 7086 is code No. X 0928.</li> <li>6. Incubate in the dark at 4 °C for 30 minutes.</li> <li>7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.</li> <li>8. Analyse on a flow cytometer.</li> </ol> <p>It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p>

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

---

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1. Davey FR, Gatter KC, Ralfkiaer E, Pulford KAF, Krissansen GW, Mason DY. Immunophenotyping of non-Hodgkin's lymphomas using a panel of antibodies on paraffin-embedded tissues. *Am J Pathol* 1987;129: 54-63.
2. Sewell WA, Cooley MA, Hegen M. NL6. CD45 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan.* New York, London: Garland Publishing Inc.; 1997. p. 499-502.
3. Akbar AN, Salmon M, Janossy G. The synergy between naive and memory T cells during activation. *Immunology Today* 1991;12:184-6.
4. Pulford KAF, Falini B, Heryet A, Gatter KC, Mason DY. N1.9. 4KB5, a new monoclonal anti-B-cell antibody for the routine diagnosis of lymphoid tissue biopsies. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. *Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England.* Oxford, New York, Tokyo: Oxford University Press; 1987. p. 828.

5. Schmidt RE. Non-lineage/natural killer section report: new and previously defined clusters. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 517-42.
6. Morimoto C. T18. CD45 cluster report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 386-9.

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Agilent Technologies Singapore (International) Pte Ltd.