Monoclonal Mouse
Anti-Human CD56/RPE
 Clone C5.9
Code No./ Code/ Code-Nr. R 7251
 Edition/ Ausgabe 01.01.03

ENGLISH	
Intended use	For in vitro diagnostic use. R 7251 is intended for use in flow cytometry. CD56 is the prototypic marker of human natural killer (NK) cells and it is also present on a subset of CD4-positive and CD8-positive T cells in peripheral blood (1). Interpretation of staining results with antibodies to CD56 must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
Synonym for antigen	Neural cell adhesion molecule (N-CAM), Leu-19, NKH1 (1).
Introduction	CD56 is a membrane glycoprotein which has multiple isoforms generated by alternative splicing from a single gene located on chromosome 11 (2). The core polypeptide is the 140 kDa isoform which can be variably glycosylated to produce mature species. The predominant isoform in NK and T cells is the transmembrane-anchored 140 kDa protein. In tissues, CD56 antigen is found in the cerebellum and cortex of the brain, and at neuromuscular junctions. This molecule is identical to the 140 kDa isoform of the neuronal cell adhesion molecule (N-CAM) (1, 3). Additionally, CD56 is also present in certain large granular lymphocyte leukaemias, small cell lung carcinomas, neural-derived tumours, myelomas, and myeloid leukaemias (1).
Reagent provided	R 7251, is a purified monoclonal mouse antibody conjugated with R-phycoerythrin (RPE). The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG2b, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.
Immunogen	Human NK cells.
Specificity	When using FACS analysis on buffy coat cells of peripheral blood as well as on gated lymphocyte populations and NK-cell cultures, Anti-CD56, C5.9, shows a positive reaction with NK-cells. The staining pattern parallels that observed with Anti-CD56, MOC-1, which has been clustered as an antibody to CD56 (2). The epitope recognized by Anti-CD56, C5.9, is heat and protease sensitive. Note, R 7251 labels CD56-positive cells with higher fluorescence intensity than Anti-CD56/RPE, Clone MOC-1, and provides a better separation between positive and negative cells.
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Staining procedure	1. Collect venous blood into a test tube containing an anticoagulant. 2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. 3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. 4. In one test tube, mix 100 µL cell suspension with 10 µL R 7251. 5. In another test tube, mix 100 µL cell suspension with an appropriate volume, e.g. 10 µL, negative control reagent. The recommended negative control for R 7251 is DakoCytomation code No. X 0951. 6. Incubate in the dark at 4 °C for 30 minutes. 7. Wash twice with PBS containing 2% BSA. Resuspend the cells in each test tube in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. 8. Analyse on a flow cytometer. It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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