Monoclonal Mouse Anti-Human CD117, c-kit, Clone 104D2

Code No./ Code/ Code-Nr. **C 7244** **APC-Conjugated**
Code No./ Code/ Code-Nr. **R 7145** **RPE-Conjugated**

ENGLISH

Intended use For in vitro diagnostic use.
C 7244 and R 7145 are intended for use in flow cytometry. Anti-CD117, 104D2, is a useful tool for the flow cytometric identification of acute myeloid leukaemia (1), and, generally, antibodies to CD117 are valuable for the classification of leukaemias as they label most myeloid, but rarely lymphoid blast cells (2, 3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Synonyms for antigen Stem cell factor receptor, mast cell growth factor receptor.

Introduction CD117 is a 145 kDa transmembrane protein which belongs to the class III receptor tyrosine kinase family. The extracellular region consists of five immunoglobulin-like domains where the second and third loops are thought to be involved in ligand binding. The natural ligand for CD117 has been termed stem cell factor, steel factor, or mast cell growth factor (3, 4). CD117 was originally identified as an acute myeloid leukaemia-associated antigen, and was also shown to be a marker for tissue mast cells, haematopoietic stem cells, and progenitor cells in normal human bone marrow. The majority of CD117-positive marrow cells (50-70%) coexpress CD34 and comprise progenitor cells and their precursors of all haematopoietic lineages. Apart from mast cells, CD117 is essentially absent from mature haematopoietic cells. Although CD117 is primarily a marker for non-lymphoid progenitor cells in the haematopoietic system, it is also found on early lymphoid progenitors, and has been reported on some pediatric T-acute lymphoblastic leukaemia (ALL) cell lines, although it seems uncommon in primary ALL. CD117 is also expressed by non-haematopoietic cells/tissues, including melanocytes, embryonic brain, and the reproductive system (3).

Reagent provided The Anti-CD117 conjugates, C 7244 and R 7145, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10⁶ Hel.92.1.7 cells).
Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
R 7145	RPE (R-Phycoerythrin)	X 0928
C 7244	APC (Allophycocyanin)	X 0968

Specificity Anti-CD117, 104D2, was included in the Sixth International Workshop and Conference on Human Leucocyte differentiation Antigens, and studies by different laboratories confirmed its reactivity with CD117 (3).

Precautions

- For professional users.
- This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
- As with any product derived from biological sources, proper handling procedures should be used.

Storage Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure

- Collect venous blood into a test tube containing an anticoagulant.
- Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
- Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
- Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD117.
- Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
- Incubate in the dark at 4 °C for 30 minutes.
- Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
- Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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

- Bahia DMM, Yamamoto M, de Lourdes M, Chauffaille LF, Kimura EYS, Bordin JO, et al. Aberrant phenotypes in acute myeloid leukaemia: a high frequency and clinical significance. *Haematologica* 2001;86:801-6.
- Bühring H-J, Ullrich A, Schaudt K, Müller CA, Busch FW. The product of the proto-oncogene c-kit (P145^{c-kit}) is a human bone marrow surface antigen of hemopoietic precursor cells which is expressed on a subset of acute non-lymphoblastic leukemic cells. *Leukemia* 1991;5:854-60.
- Ashman LK, Cambareri AC, Nguyen L, Bühring H-J. CR5. CD117 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference*; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 816-8.
- Bühring H-J, Ashman LK, Gattei V, Kniep B, Larregina A, Pinto A, et al. CR2.7. Stem-cell factor receptor (p145(c-kit)) summary report (CD117). In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference*; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 1882-8.

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