Monoclonal Mouse Anti-Human CD79 α cy/APC, Clone HM57
Monoclonal Mouse Anti-Human CD79 α cy/RPE, Clone HM57

Code C7252
Code R7159

ENGLISH

Intended use	For in vitro diagnostic use. C7252 and R7159 are intended for use in flow cytometry. CD79 α is considered as a pan-B cell antigen, widely used in the analysis of B cell malignancies (5). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.									
Synonyms for antigen	Ig α , IgM- α , MB1 (4, 6)									
Summary and explanation	<p>The CD79α polypeptide forms, together with the CD79β polypeptide, a disulphide-linked, transmembrane heterodimer, which is non-covalently associated with membrane-bound immunoglobulins on B cells (1, 4). These 3 components constitute the B cell antigen receptor complex. Antigen binding by the immunoglobulin in this complex leads to CD79α/CD79β-mediated signal transduction resulting in activation of multiple biochemical pathways (4). CD79α/CD79β dimer may thus be considered to be the B-cell equivalent to CD3 on T cells (1). The CD79α and CD79β have a molecular mass of 32-33 kDa and 37-39 kDa, respectively, and their cytoplasmic determinants are named CD79αcy and CD79βcy (2).</p> <p>In Normal cells CD79α expression are restricted B cell lineage. CD79αcy is found in the cytoplasm of pro B- cells before Ig heavy-chain rearrangement, and is extinguished in terminally differentiated plasma cells. CD79α is expressed in hairy cell leukaemia and in the majority of all low-grade B-cell leukaemias and lymphomas, such as lymphoplasmacytic, follicular, mantle cell, marginal zone and Burkitt's lymphomas. Generally, the CD79α antigen is not found in acute myeloid leukaemia, but the FAB M3 class is a significant exception with more than 50% of cases displaying the antigen in the cytoplasm of malignant cells. Approximately 45% of T- cell ALL cases co-express CD3 and CD79α, and this co-expression phenotype may be clinically relevant (4).</p>									
Reagent provided	<p>The Anti-CD79αcy, conjugates, C 7252 and R 7159, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 μL of conjugate for up to 10⁹ leucocytes from normal human peripheral blood).</p> <p><u>Isotype</u>: IgG1, kappa. <u>Conjugate concentration mg/L</u>: See label on vial.</p> <table border="1" data-bbox="381 943 1240 1068"> <thead> <tr> <th>Antibody Code</th> <th>Fluorochrome</th> <th>Negative Control Code</th> </tr> </thead> <tbody> <tr> <td>R7159</td> <td>RPE (R-Phycoerythrin)</td> <td>X0928</td> </tr> <tr> <td>C7252</td> <td>APC (Allophycocyanin)</td> <td>X0968</td> </tr> </tbody> </table>	Antibody Code	Fluorochrome	Negative Control Code	R7159	RPE (R-Phycoerythrin)	X0928	C7252	APC (Allophycocyanin)	X0968
Antibody Code	Fluorochrome	Negative Control Code								
R7159	RPE (R-Phycoerythrin)	X0928								
C7252	APC (Allophycocyanin)	X0968								
Immunogen	Synthetic peptide sequence comprising amino acids 202-216 of the CD79 α protein. This oligopeptide represents an intra-cytoplasmic C-terminal part of the protein (1, 3).									
Specificity	<p>The HM57 antibody was clustered as Anti-CD79αcy at the Fifth International Workshop on Human Leucocyte Differentiation Antigens in Boston in 1993 (2). In tissue sections Anti-CD79αcy, HM57, labels B-cells from several mammalian species, i.e. calf, guinea pig, horse, monkey, mouse, rat and swine (1).</p> <p>In normal tissue Anti-CD79αcy, HM57, detects CD79α in cells restricted to B-cell lineage. In hairy cell leukaemia and in the majority of all low-grade B cell leukaemias and lymphomas, such as lymphoplasmacytic, follicular, mantle cell, marginal zone and Burkitt's lymphomas, CD79αcy can be detected using Anti-CD79αcy, HM57 (4).</p>									
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations. 									
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
Staining procedure	<ol style="list-style-type: none"> Transfer 50 μL (up to 10⁶ cells) of the cell suspension to be analysed (whole blood, bone marrow or mononuclear cells) to a test tube. Add 100 μL Dako IntraStain, Reagent A (Fixation), code No. K2311. Mix gently by using a vortex mixer to ensure that the cells are in suspension. Incubate at room temperature for 15 minutes. Add 2 mL PBS and mix gently by using a vortex mixer. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 μL of fluid. Mix thoroughly by using a vortex mixer to ensure that the cells are in suspension and add 100 μL Dako IntraStain, Reagent B (Permeabilization), code No. K2311. Add 10 μL fluorochrome-conjugated Anti-CD79αcy. Mix gently by using a vortex mixer to ensure that the cells are in suspension. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at room temperature for 15 minutes. Repeat steps 4 and 5. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. Analyse on a flow cytometer 									

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
1. Mason DY, Cordell JL, Tse AGD, van Dongen JJ M, van Noesel CJM, Micklem K. The IgM-associated protein mB 1 as a marker of normal and neoplastic B cells. *J Immunol* 1991;147:2474-82.
2. Engel P, Wagner N, Tedder TF. B23. CD79 workshop report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 1. p. 667-70.*
3. Cordell JL, Jones M, Brown M, Gatter KC, van Dongen JJM, Mason DY. B23.4. Studies of the mb1 polypeptide - a major new B cell marker. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 1. p. 676-7.*
4. Chu PG, Arber DA. CD79: A review. *Appl Immunohistochem Mol Morphology* 2001;9:97-106.
5. Verschuren MCM, Comans-Bitter WM, Mason DY, Drexler HG, van Dongen JJM. B23.4. Expression of mb1 and B29 proteins in human haemopoietic malignancies and cell lines. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 1. p. 677-79.*
6. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. *Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 828.*

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Agilent Technologies Singapore (International) Pte Ltd.