Monoclonal Mouse Anti-Human Myeloperoxidase/APC, Clone MPO-7
Monoclonal Mouse Anti-Human Myeloperoxidase/FITC, Clone MPO-7
Monoclonal Mouse Anti-Human Myeloperoxidase/RPE, Clone MPO-7
Monoclonal Mouse Anti-Human Myeloperoxidase/PerCP-Cy5.5, Clone MPO-7

Code C7246
Code F0714
Code R7209
Code PR704

ENGLISH

Intended use	For in vitro diagnostic use. C7246, F0714, R7209 and PR704 are intended for use in flow cytometry. The antibody labels myeloperoxidase-expressing cells in normal and neoplastic tissue (4). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
Summary and explanation	Human myeloperoxidase (MPO) is a 150 kDa tetrameric polypeptide consisting of 2 covalently linked subunits of approximately 59 kDa and 2 subunits of approximately 13 kDa. MPO is a lysosomal enzyme present abundantly in the azurophilic granules of neutrophils and at lower concentration in the granules of monocytes. MPO is involved in the microbicidal action of these cells (1).
Reagent provided	The Anti-MPO conjugates, C7246, F0714, R7209 and PR704, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains conjugate for 100 tests of up to 10 ⁶ leucocytes from normal human peripheral blood. See table 1 for volume of conjugate per test. <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration:</u> See label on vial.

Table 1

Antibody Dako Code	Fluorochrome	Volume of conjugate per test	Isotype Reagent Dako Code
F0714	FITC (Fluorescein Isothiocyanate Isomer 1)	10 µL	X0927
R7209	RPE (R-Phycoerythrin)	10 µL	X0928
C7246	APC (Allophycocyanin)	10 µL	X0968
PR704	PerCP-Cy5.5 (Peridinin Chlorophyll Protein Cyanine 5.5)	5 µL	No Dako Code*

*It is recommended to use a non-reactive monoclonal antibody of the same isotype and conjugated to PerCP-Cy5.5.

Immunogen	Human MPO isolated from CML cells (2).
Specificity	Anti-MPO, MPO-7, strongly labels the cytoplasm of neutrophil granulocytes in peripheral blood. Monocytes are weakly positive while lymphocytes are unreactive. In bone marrow, cells of the granulocyte series are labelled. Anti-MPO, MPO-7, also labels the human promyelocyte cell line HL-60, while other myeloid cell lines such as K-562 and KG-1 are consistently unreactive (3). In acute myeloid leukemia, anti-MPO, MPO-7, is a useful reagent (4) since it labels a large proportion of cases (3, 5).
Precautions	<ol style="list-style-type: none"> For in vitro diagnostic use. For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Staining procedure	<ol style="list-style-type: none"> Transfer 50 µL (up to 10⁶ cells) of the cell suspension to be analysed (whole blood, bone marrow or mononuclear cells) to a test tube. Add 100 µL Dako IntraStain Reagent A (Fixation), Code K2311. Mix gently by using a vortex mixer to ensure that the cells are in suspension. Incubate at room temperature for 15 minutes. Add 2 mL PBS and mix gently by using a vortex mixer. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 µL of fluid. Mix thoroughly by using a vortex mixer to ensure that the cells are in suspension and add 100 µL Dako IntraStain Reagent B (Permeabilization), Code K2311. Add the volume of anti-MPO conjugate stated in table 1. Mix gently by using a vortex mixer to ensure that the cells are in suspension.

7. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
8. Incubate in the dark at room temperature for 15 minutes.
9. Repeat steps 4 and 5.
10. Resuspend the pellet in an appropriate fluid for flow cytometric analysis, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
11. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
1. Klebanoff SJ. Myeloperoxidase. Proc Assoc Am Physicians 1999;111:383-89.
2. Olsson I, Olofsson T, Odeberg H. Myeloperoxidase-mediated iodination in granulocytes. Scand J Haemat 1972;9:483-91.
3. Nakase K, Sartor M, Bradstock K. Detection of myeloperoxidase by flow cytometry in acute leukemia. Cytometry 1998;34:198-202.
4. Kappelmayer J, Gratama JW, Karázi É, Menéndez P, Ciudad J, Rivas R, et al. Flow cytometric detection of intracellular myeloperoxidase, CD3 and CD79a. Interaction between monoclonal antibody clones, fluorochromes and sample preparation protocols. J Immunol Methods 2000;242:53-65.
5. Lanza F, Latorraca A, Moretti S, Castagnari B, Ferrari L, Castoldi G. Comparative analysis of different permeabilization methods for the flow cytometry measurement of cytoplasmic myeloperoxidase and lysozyme in normal and leukemic cells. Cytometry 1997;30:134-44.

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Agilent Technologies Singapore (International) Pte Ltd.