Monoclonal Mouse Anti-Human CD64, Fc Gamma Receptor I/APC, Clone 10.1
Monoclonal Mouse Anti-Human CD64, Fc Gamma Receptor I/RPE, Clone 10.1
Monoclonal Mouse Anti-Human CD64, Fc Gamma Receptor I/RPE-Cy5, Clone 10.1

Code C7278
Code R7219
Code C7220

ENGLISH													
Intended use	For in vitro diagnostic use. C7278, R7219 and C7220 are intended for use in flow cytometry. Anti-CD64, 10.1, is intended for use in identifying cells expressing CD64. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.												
Synonyms for antigen	FC γ RI, FCR I (1).												
Introduction	CD64 is the high-affinity IgG Fc receptor belonging to the immunoglobulin supergene family (1). CD64 is a 72 kDa glycoprotein which binds monomeric IgG (human IgG1>IgG3>>IgG4>>>IgG2 and mouse IgG3>IgG2a>>>IgG1, IgG2b). Three genes: <i>FcγR1A</i> , <i>FcγR1B</i> , and <i>FcγR1C</i> , which have >98 percent identity at the nucleotide level, have been identified and mapped to chromosome 1, band q21.1 (2). The <i>FcγR1A</i> gene, which specifies the high-affinity IgG receptor, encodes three C2-type Ig-like extracellular domains, a 21 amino acid transmembrane region, and a charged cytoplasmic tail of 61 amino acids (2, 3). CD64 is constitutively expressed on monocytes, macrophages and blood dendritic cells, and at very low levels on neutrophils. It can be induced to higher expression on neutrophils and eosinophils by interferon (IFN) γ and granulocyte colony-stimulating factor (G-CSF) (2). CD64 is expressed in acute myeloid leukaemia with monocytic differentiation (M4 and M5) and in acute promyelocytic leukaemia (M3) (4, 5). CD64 mediates endocytosis, phagocytosis, antibody-dependent cellular toxicity, cytokine release, and superoxide generation (2).												
Reagent provided	The Anti-CD64 conjugates, C7278, R7219 and C7220, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 μ L of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.												
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Control Reagent Code No.</th> </tr> </thead> <tbody> <tr> <td>C7278</td> <td>APC (Allophycocyanin)</td> <td>X0968</td> </tr> <tr> <td>R7219</td> <td>RPE (R-Phycoerythrin)</td> <td>X0928</td> </tr> <tr> <td>C7220</td> <td>RPE-Cy5 (R-Phycoerythrin-Cyanine 5)</td> <td>X0955</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Control Reagent Code No.	C7278	APC (Allophycocyanin)	X0968	R7219	RPE (R-Phycoerythrin)	X0928	C7220	RPE-Cy5 (R-Phycoerythrin-Cyanine 5)	X0955
Antibody Code No.	Fluorochrome	Control Reagent Code No.											
C7278	APC (Allophycocyanin)	X0968											
R7219	RPE (R-Phycoerythrin)	X0928											
C7220	RPE-Cy5 (R-Phycoerythrin-Cyanine 5)	X0955											
Immunogen	Human rheumatoid synovial fluid cells and monocytes (6).												
Specificity	Anti-CD64, 10.1, was included in the Fifth and Sixth International Workshops and Conferences on Human Leucocyte Differentiation Antigens. Studies by a number of laboratories confirmed its reactivity with the CD64 antigen (2, 7). Anti-CD64, 10.1, reacts with monocytes and the myeloid cell lines HL-60, U937, THP-1, and Mono Mac 6. It shows no reaction with resting and activated B cells, T cells, NK cells, eosinophils, platelets, fibroblasts, endothelial cells, and epithelial cells (2, 6). Anti-CD64, 10.1, reacts with an epitope on Fc γ R1 near to the binding site for the Fc region of IgG and blocks the binding of human IgG, rabbit IgG and mouse IgG2a to Fc γ R1 (6).												
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.												
Storage	Store in the dark at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.												
Staining procedure	<ol style="list-style-type: none"> Transfer 100 μL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube. Add 10 μL of fluorochrome-conjugated Anti-CD64 and mix gently with a vortex mixer. The 10 μL is a guideline only; the optimal volume should be determined by the individual laboratory. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes or at room temperature (20–25 °C) for 15–30 minutes. Add 100 μL of DakoCytomation Uti-Lyse™ (code Nos. S3325 or S3350) Reagent A to each sample and mix gently with a vortex mixer. Incubate for 10 minutes at room temperature in the dark. Add 1 mL of DakoCytomation Uti-Lyse™ Reagent B to each sample and mix gently with a vortex mixer. Incubate for 10 minutes at room temperature in the dark. If another lysing reagent is used in steps 5 and 6, please follow the recommendations for that reagent. Centrifuge at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 μL of fluid. Add 2 mL 0.01 mol/L PBS containing 2% bovine serum albumin and resuspend the cells by using a vortex mixer. Repeat step 7. Resuspend pellet in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS. The PBS should contain 1% paraformaldehyde (fixative) if samples are not analysed the same day. Analyse on a flow cytometer or store at 2–8 °C in the dark until analysis. Samples can be run up to 24 hours after lysis. 												

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Abnormal numbers of cells expressing the target antigens or aberrant expression levels of the antigens can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern for the antigens and their relationship to expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour analysis for the comprehensive analysis of neoplastic cells of leukaemia and lymphoma.

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
References/ Références/ Literatur	<ol style="list-style-type: none"> Guyre P, Sun W. CD Guide. CD64. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 815. Sun W, O'Shea JK, Guyre PM. MC9. CD64 Workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 988-90. van de Winkel JGJ, Capel PJA. Human IgG Fc receptor heterogeneity: molecular aspects and clinical implications [review]. Immunol Today 1993;14:215-21. Krasinskas AM, Wasik MA, Kamoun M, Schretzenmair R, Moore J, Salhany KE. The usefulness of CD64, other monocyte-associated antigens, and CD45 gating in the subclassification of acute myeloid leukemias with monocytic differentiation. Am J Clin Pathol 1998;110:797-805. Repp R, Schaekel U, Helm G, Thiede C, Soucek S, Pascheberg U, et al. Immunophenotyping is an independent factor for risk stratification in AML. Cytometry 2003;53B:11-9. Dougherty GJ, Selvendran Y, Murdoch S, Palmer DG, Hogg N. The human mononuclear phagocyte high-affinity Fc receptor, FcR1, defined by a monoclonal antibody, 10.1. Eur J Immunol 1987;17:1453-9. Guyre PM, Von dem Borne AEG. M12. CD64 cluster workshop report. In: Schlossman SF, Bousmell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 1. p. 874-5. van Vugt MJ, van den Herik-Oudijk IE, van de Winkle JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). Blood 1996;88:2358-61.
Explanation of symbols/ Légende des symboles/ Erläuterung der Symbole	

REF	Catalogue number Référence du catalogue Bestellnummer	 2°C–8°C	Temperature limitation Limites de température Zulässiger Temperaturbereich		Use by Utiliser jusque Verwendbar bis
IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum		Keep away from sunlight (consult storage section) Conserver à l'écart du soleil (se reporter à la section conservation) Lichtgeschützt lagern (siehe Abschnitt zur Lagerung)		Manufacturer Fabricant Hersteller
	Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	LOT	Batch code Code du Lot Chargenbezeichnung		