MultiMix™
Triple-Colour Reagent
Anti-Human CD3/FITC
Anti-Human CD4/RPE
Anti-Human CD45/RPE-Cy5
Code No./ Code/ Code-Nr. TC646

ENGLISH	
Intended use	For in vitro diagnostic use. TC646 is intended for use in flow cytometry. CD3 is a pan-T-cell restricted antigen which is a valuable marker for normal and neoplastic T cells. Measurement of the relative as well as the absolute number of CD4-positive cells in peripheral blood is critically important in assessing the immunological status of a patient. In flow cytometry, anti-CD45, together with a panel of other antibodies, is considered essential for the initial evaluation of chronic lymphoproliferative disorders and acute leukaemias (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
Reagent provided	TC646 comprises the following three, carefully matched, fluorescent antibodies: Monoclonal Mouse Anti-Human CD3, Clone UCHT1, conjugated with fluorescein isothiocyanate isomer 1 (FITC). Monoclonal Mouse Anti-Human CD4, Clone MT310, conjugated with R-phycoerythrin (RPE). Monoclonal Mouse Anti-Human CD45, Clone T29/33, conjugated with R-phycoerythrin-Cyanine 5 (RPE-Cy5). The three conjugates in TC646 have been produced from purified monoclonal mouse antibodies of isotype IgG1, kappa. TC646 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 50 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood).
Specificity	Anti-CD3, UCHT1, was included in the First and Third International Workshops and Conferences on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD3 (2). Anti-CD3, UCHT1, reacts with the ε-chain of CD3 (3). Anti-CD4, MT310, was included in the Second International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD4 (4). Anti-CD45, T29/33, was included in the Third and Fourth International Workshops and Conferences on Human Leucocyte Differentiation Antigens. At the Third Workshop, the antibody was clustered as a CD45 antibody reacting with all known isotypes of the CD45 family, also called the leucocyte common antigen family (5). At the Fourth Workshop, expression of CD45 in a range of haematopoietic cell lines, and lack of CD45 in non-haematopoietic cell lines, were demonstrated. Notably, CD45 was expressed at a very low level in myeloma cells and in the U-266 plasmacytoid cell line (6).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Staining procedure	1. Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube. 2. Add 10 µL of TC646 and mix gently by using a vortex mixer. 3. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes. 4. Add 100 µL of Uti-Lyse™ Reagent A (DakoCytomation code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark. 5. Add 1 mL of Uti-Lyse™ Reagent B (DakoCytomation code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark. 6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube. 7. Add 2 mL of PBS (DakoCytomation code No. S3024) to the tube and resuspend the cells by using a vortex mixer. 8. Repeat step 6. 9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS. 10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.
Procedural notes	Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control. Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory. It is optional to include a control reagent test tube. The control reagent should match the isotypes and fluorochromes of the conjugated antibody reagent. The recommended Triple-Colour Control Reagent for TC646 is DakoCytomation code No. X0964.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as DakoCytomation EasyLyse™, code No. S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis.

Product-specific limitations It has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (7).

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
- Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematology neoplasias: results of an international consensus meeting. *Cytometry* 2001;46:23-7.
- McMichael AJ, Gotch FM. T-cell antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. *Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference*; 1986 Sept 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 31-62.
- Tunnacliffe A, Olsson C, Traunecker A, Krissansen GW, Karjalainen K, de la Hera A. T3.2. The majority of CD3 epitopes are conferred by the epsilon chain. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. *Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference*; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 295-6.
- Miedema FT, FG. Melief, CJM. T cell-dependent immunoglobulin synthesis in the human system. Studies with T cell-specific monoclonal antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. *Leukocyte typing II. Proceedings of the 2nd International Workshop on Human Leukocyte Differentiation Antigens*; 1984 Sept 17-20; Boston, USA. New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume I. p. 213-22.
- Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1, and leucocyte common antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. *Leukocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference*; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 788-803.
- Schwinger R. N7. Cluster report: CD45/CD45R. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. *Leukocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference*; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 628-34.
- van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). *Blood* 1996;88:2358-61.

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