Isotype Reagent	
Mouse IgG1/APC	Code X0968
Mouse IgG1/FITC	Code X0927
Mouse IgG1/PerCP	Code X7909
Mouse IgG1/PerCP-Cy5.5	Code X7914
Mouse IgG1/RPE	Code X0928
Mouse IgG1/RPE-Cy5	Code X0955

ENGLISH

Intended use	For in vitro diagnostic use. X0968, X0927, X7909, X0928, X0955 and X7914 are intended for use in flow cytometry. These reagents have been produced especially as negative controls for Dako fluorochrome-conjugated monoclonal mouse antibodies of isotype IgG1.
Summary and explanation	In flow cytometry, isotype-specific negative control reagents typically are fluorochrome-conjugated monoclonal antibodies directed against non-human antigens. These reagents measure the amount of conjugate binding that is unrelated to the target antigen of the human cell. This non-immunological binding may be caused by Fc receptors on the cells or by the physico-chemical characteristics of the fluorochrome conjugate. A conjugated monoclonal mouse antibody selected as an isotypic control should be of the same isotype as the conjugated monoclonal mouse antibody used as the test reagent and should be conjugated with the same fluorochrome.
Reagent provided	X0968, X0927, X7909, X0928, X0955 and X7914 are all produced from a purified monoclonal mouse antibody clone DAK-GO1, conjugated with fluorescein isothiocyanate 1 (FITC), R-phycoerythrin (RPE), peridinin chlorophyll protein (PerCP), R-phycoerythrin-Cyanine5 (RPE-Cy5), peridinin chlorophyll protein-Cyanine5.5 (PerCP-Cy5.5) or cross-linked allophycocyanin (APC), respectively. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.
Specificity	The antibody specificity of the mouse IgG1 produced by clone DAK-GO1 is directed towards <i>Aspergillus niger</i> glucose oxidase, an enzyme which is neither present nor inducible in mammalian tissues.
Precautions	<ol style="list-style-type: none"> 1. For professional users. 2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Dilution guideline	In flow cytometry, the conjugated negative controls may be used at a volume of 10 µL per test when adjusted to the same conjugate concentration as the test antibody. This is a guideline only. The optimal amount of the negative control conjugate may vary, e.g. depending on the specimen and preparation method, and should be determined by each individual laboratory. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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Agilent Technologies Singapore (International) Pte Ltd.