Ig Reagent
Rabbit F(ab')₂/APC **Code X0998**
Rabbit F(ab')₂/FITC **Code X0929**
Rabbit F(ab')₂/RPE **Code X0930**

ENGLISH

Intended use	For in vitro diagnostic use. X0998, X0929 and X0930 are intended for use in flow cytometry. These reagents have been produced especially as negative controls for Dako APC, FITC-, or RPE-conjugated F(ab') ₂ fragments of rabbit antibodies.
Summary and explanation	Negative control reagents for F(ab') ₂ fragments of fluorochrome-conjugated rabbit antibodies measure the amount of antibody binding that is unrelated to the target antigen present on the human cell. This non-immunological binding may be caused by the physico-chemical characteristics of the fluorochrome conjugate.
Reagent provided	X0998, X0929 and X0930 are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Conjugate concentration g/L:</u> See label on vial.
Preparation	<ol style="list-style-type: none"> 1. The immunoglobulin (Ig) fraction is isolated from serum of healthy, non-immunized rabbits (Danish Whites) by salting out and ion exchange chromatography. 2. The Ig fraction is solid-phase absorbed with human plasma. This reduces the non-specific background, and makes the reagent a better control for solid-phase absorbed antibodies. 3. The immunoglobulin molecules are then digested by pepsin and the F(ab')₂ fraction isolated by gel filtration. 4. Finally, the F(ab')₂ fraction is conjugated with allophycocyanin (APC), fluorescein isothiocyanate isomer 1 (FITC), or R-phycoerythrin (RPE).
Precautions	<ol style="list-style-type: none"> 1. For professional users. 2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.
Application	In flow cytometry, X0998, X0929 and X0930 may be used at a volume of 10 µL per test when adjusted to the same conjugate concentration as the test antibody. This is a guideline only. The optimal amount of the negative control may vary, e.g. depending on the specimen and preparation method and should be determined by each individual laboratory. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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Agilent Technologies Singapore (International) Pte Ltd.