**MultiMix™**  
**Dual-Colour Control Reagent**  
**Rabbit F(ab')<sub>2</sub>/FITC**  
**Rabbit F(ab')<sub>2</sub>/RPE**  
**Code No./ Code/ Code-Nr. X0935**

ENGLISH	
<b>Intended use</b>	For in vitro diagnostic use. X0935 is intended for use in flow cytometry. X0935 has been produced especially as a flow cytometry control reagent for DakoCytomation dual-colour FITC- and RPE-conjugated F(ab') <sub>2</sub> fragments of rabbit antibodies.
<b>Reagent provided</b>	X0935 comprises the following two, carefully matched, fluorescent rabbit F(ab') <sub>2</sub> reagents: F(ab') <sub>2</sub> fragment of rabbit antibodies, conjugated with fluorescein isothiocyanate isomer 1 (FITC). F(ab') <sub>2</sub> fragment of rabbit antibodies, conjugated with R-phycoerythrin (RPE). X0935 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 50 tests (10 µL of conjugate for up to 10 <sup>6</sup> leucocytes from normal human peripheral blood).
<b>Preparation</b>	<ol style="list-style-type: none"> <li>The immunoglobulin (Ig) fraction is isolated from serum of healthy, non-immunized rabbits (Danish Whites) by salting out and ion exchange chromatography.</li> <li>The Ig fraction is solid-phase absorbed with human plasma. This reduces the non-specific background, and makes the reagent a better control for solid-phase absorbed antibodies.</li> <li>The immunoglobulin molecules are then degraded by pepsin and the F(ab')<sub>2</sub> fraction isolated by gel filtration.</li> <li>Finally, the F(ab')<sub>2</sub> fraction is conjugated with FITC or RPE.</li> </ol>
<b>Precautions</b>	<ol style="list-style-type: none"> <li>For professional users.</li> <li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>As with any product derived from biological sources, proper handling procedures should be used.</li> </ol>
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.
<b>Dilution guideline</b>	In flow cytometry, X0935 may be used at a volume of 10 µL per test. This is a guideline only. The optimal amount may vary, e.g. depending on the specimen and preparation method, and should be determined by each individual laboratory. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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