**MultiMix™**
**Dual-Colour Control Reagent**
**Mouse IgG1/FITC**
**Mouse IgG2a/RPE**
**Code No./ Code/ Code-Nr. X0949**
**ENGLISH**

<b>Intended use</b>	For in vitro diagnostic use. X0949 is intended for use in flow cytometry. X0949 has been produced especially as a flow cytometry control reagent for Dako dual-colour monoclonal mouse antibodies of isotype IgG1 and IgG2a.
<b>Reagent provided</b>	X0949 comprises the following two, carefully matched, fluorescent monoclonal mouse IgG1 and IgG2a reagents: Monoclonal mouse IgG1, Clone DAK-GO1, conjugated with fluorescein isothiocyanate isomer 1 (FITC). Monoclonal mouse IgG2a, Clone DAK-GO5, conjugated with R-phycoerythrin (RPE). The two conjugates in X0949 have been produced from purified monoclonal mouse antibodies. X0949 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 50 tests (10 µL of conjugate for up to 10 <sup>6</sup> leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa (Clone DAK-GO1) and IgG2a, kappa (Clone DAK-GO5).
<b>Specificity</b>	The antibody specificity of mouse IgG1 and IgG2a produced by clones DAK-GO1 and DAK-GO5, respectively, is directed towards <i>Aspergillus niger</i> glucose oxidase, an enzyme which is neither present nor inducible in mammalian tissues.
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.
<b>Dilution guideline</b>	In flow cytometry, X0949 may be used at a volume of 10 µL. This is a guideline only. The optimal amount may vary, e.g. depending on the specimen and preparation method, and should be determined by each individual laboratory. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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Agilent Technologies Singapore (International) Pte Ltd.