**MultiMix™**  
**Triple-Colour Control Reagent**  
**Mouse IgG1/FITC**  
**Mouse IgG1/RPE**  
**Mouse IgG1/RPE-Cy5**  
**Code No./ Code/ Code-Nr. X0956**

ENGLISH	
<b>Intended use</b>	For in vitro diagnostic use. X0956 is intended for use in flow cytometry. X0956 has been produced especially as a flow cytometry control reagent for DakoCytomation triple-colour monoclonal mouse antibodies of isotype IgG1.
<b>Reagent provided</b>	X0956 comprises the following three, carefully matched, fluorescent monoclonal mouse IgG1 reagents: Monoclonal mouse IgG1, clone DAK-GO1, conjugated with fluorescein isothiocyanate isomer 1 (FITC). Monoclonal mouse IgG1, clone DAK-GO1, conjugated with R-phycoerythrin (RPE). Monoclonal mouse IgG1, clone DAK-GO1, conjugated with R-phycoerythrin-Cyanine 5 (RPE-Cy5). The three conjugates in X0956 have been produced from purified monoclonal mouse antibodies. X0956 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 50 tests (10 µL of conjugate for up to 10 <sup>6</sup> leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa for all three conjugates in X0956.
<b>Specificity</b>	The antibody specificity of the mouse IgG1 produced by clone DAK-GO1 is directed towards <i>Aspergillus niger</i> glucose oxidase, an enzyme which is neither present nor inducible in mammalian tissues.
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.
<b>Dilution guidelines</b>	In flow cytometry, X0956 may be used at a volume of 10 µL per test. This is a guideline only. The optimal amount may vary, e.g. depending on the specimen and preparation method, and should be determined by each individual laboratory. As it has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (1), X0956 is an especially useful control reagent, since it aids in the discrimination between immunological and non-immunological reactivity of the test conjugate. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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