MultiMix™
Triple-Colour Control Reagent
Rabbit F(ab')₂/FITC
Rabbit F(ab')₂/RPE
Mouse IgG1/RPE-Cy5
Code No./ Réf./ Code-Nr. X0957

ENGLISH	
Intended use	For in vitro diagnostic use. X0957 is intended for use in flow cytometry. X0957 has been produced especially as a flow cytometry control reagent for DakoCytomation triple-colour rabbit F(ab') ₂ and monoclonal mouse antibodies of isotype IgG1.
Reagent provided	X0957 comprises the following three, carefully matched, fluorescent reagents: F(ab') ₂ fragment of rabbit antibodies, conjugated with fluorescein isothiocyanate isomer 1 (FITC). F(ab') ₂ fragment of rabbit antibodies, conjugated with R-phycoerythrin (RPE). Monoclonal mouse IgG1, Clone DAK-GO1, conjugated with R-phycoerythrin-Cyanine 5 (RPE-Cy5).
Preparation	Rabbit F(ab') ₂ /FITC and Rabbit F(ab') ₂ /RPE have been prepared as follows: <ol style="list-style-type: none"> The rabbit immunoglobulin (Ig) fraction is isolated from serum of healthy, non-immunized rabbits (Danish Whites) by salting out and ion exchange chromatography. The Ig fraction is solid-phase absorbed with human plasma. This reduces the non-specific background, and makes the reagent a better control for solid-phase absorbed antibodies. The immunoglobulin molecules are then degraded by pepsin and the F(ab')₂ fraction isolated by gel filtration. Finally, the F(ab')₂ fraction is conjugated with FITC or RPE. <p>The mouse IgG1 antibody has been produced from purified monoclonal antibody, and conjugated with RPE-Cy5. The mouse antibody produced by clone DAK-GO1 is directed towards <i>Aspergillus niger</i> glucose oxidase, an enzyme which is neither present nor inducible in mammalian tissues.</p> <p>X0957 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 50 tests (10 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).</p>
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.
Dilution guideline	In flow cytometry, X0957 may be used at a volume of 10 µL per test. This is a guideline only. The optimal amount may vary, e.g. depending on the specimen and preparation method, and should be determined by each individual laboratory. As it has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (1), X0957 is an especially useful control reagent, since it aids in the discrimination between immunological and non-immunological reactivity of the test conjugate. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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