**Monoclonal Mouse
Anti-Human CD57/FITC**
Clone TB01
Code No./ Code/ Code-Nr. F 7270

ENGLISH

Intended use For in vitro diagnostic use.
F 7270 is intended for use in flow cytometry. The antibody is intended for use in identification of cells expressing CD57. CD57 was originally identified as a natural killer (NK) cell marker, but it is not restricted to NK cells, being shared by cells belonging to different lineages (1). Antibodies to CD57 are useful for evaluation of lymphoproliferative disorders, especially large granular lymphocyte leukaemia, when used together with a panel of other antibodies (2-4). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Synonyms for antigen HNK-1, Leu-7 (1, 2).

Introduction CD57 is a carbohydrate antigen present on a variety of both haematopoietic and neuroectodermal-derived cells (3). The important carbohydrate component could account for the different molecular weights that have been attributed to CD57 when the molecule was immunoprecipitated from different cellular sources (1). In brain, the CD57 carbohydrate is associated with myelin-associated glycoprotein (MAG), NCAM and other proteins (5). In blood, the proteins associated with the antigen have not been characterized in detail.
CD57 is present on subsets of NK cells, CD8-positive lymphocytes, and on a small percentage of CD4-positive/CD45R0-positive T cells (1, 2). In contrast, CD57 is not present on red blood cells, granulocytes, monocytes, platelets, cord blood lymphocytes, immature resting or activated thymocytes (1, 5). In neoplastic tissues, CD57 is present in large granular lymphocyte leukaemia to various degrees (4). Further, CD57 is present in Ewing's sarcomas, melanomas, neuroblastomas, neuroectodermal tumours, neurofibromas and small-cell lung carcinomas (2, 3).

Reagent provided The Anti-CD57 conjugate, F 7270, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).
Isotype: IgM, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 7270	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0934

Immunogen A pool of human neuroblastoma cell lines (6).

Specificity/reactivity Anti-CD57, TB01, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD57 (1).
The antibody recognizes the CD57 carbohydrate epitope present on several surface molecules. In the nervous system this epitope is particularly prominent on myelin-associated glycoprotein (MAG) (1). The epitope for TB01 is neuraminidase- and bromelaine-sensitive, but trypsin-resistant (1).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. This product contains material of animal origin, therefore the product should be handled as potentially infectious.

Storage Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure

1. Collect venous blood into a test tube containing an anticoagulant.
2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
4. Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD57.
5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
6. Incubate in the dark at 4 °C for 30 minutes.
7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.

8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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

1. Funaro A, Malavasi F. NK5. CD57 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 274-6.
2. Poggi A, CD guide In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 1156.
3. Funaro A, Cafforio P, Horenstein A, Magrini E, Ardeleanu C, Guena et al. Human CD57, a link molecule between leucocytes and neural cells. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 2. p. 1435-6.
4. Macey MG, McCarthy DA, Newland AC. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 276-7.
5. Schubert J, Lanier LL, Schmidt RE. N17 Cluster report: CD57. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 711-14.
6. Funaro A, Cafforio P, Horenstein A, Magrini E, Ardeleanu C, Guena et al. NK3.6 Epitope analysis of human CD57 by means of a panel of newly generated high affinity murine mAb. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 2. p. 1412-14.

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Agilent Technologies Singapore (International) Pte Ltd.