**MultiMix™**  
**Triple-Colour Reagent**  
**Anti-Human CD103/FITC**  
**Anti-Human CD11c/RPE**  
**Anti-Human CD19/APC**  
**Code TC665**

**ENGLISH**

**Intended use**

For in vitro diagnostic use.

TC665 is intended for use in flow cytometry. The reagent is intended for use in identification of cells expressing CD103, CD11c and CD19.

CD103 is expressed on mucosal T cells and on activated CD8-positive T cells. In several malignant conditions, such as T-cell lymphomas and hairy cell leukaemias, the cells express CD103. Antibodies to CD103 are valuable for the evaluation of chronic B-cell leukaemias and T-cell lymphomas (1, 2).

CD11c is expressed on a variety of cells including granulocytes, monocytes, macrophages, natural killer (NK)-cells, dendritic cells and neoplastic cells. Antibodies to CD11c are useful for the initial evaluation of B-cell lymphoproliferative disorders, e.g. hairy cell leukaemia and B-cell chronic lymphocytic leukaemia (1, 3, 4).

CD19 is the broadest lineage-specific surface marker for B cells and it is present on the surface of virtually all B lymphocytes, including early B cell progenitors (5). CD19 expression is maintained in B-lineage cells that have undergone neoplastic transformation (6). Antibodies to CD19 are considered essential for the initial evaluation of acute and chronic lymphoproliferative disorders (7).

Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Reagent provided**

TC665 comprises the following three matched fluorescent antibodies:

Monoclonal Mouse Anti-Human CD103, Clone Ber-ACT8, conjugated with fluorescein isothiocyanate isomer 1 (FITC).

Monoclonal Mouse Anti-Human CD11c, Clone KB90, conjugated with R-phycoerythrin (RPE).

Monoclonal Mouse Anti-Human CD19, Clone HD37, conjugated with allophycocyanin (APC).

The three conjugates in TC665 have been produced from purified monoclonal mouse antibodies of isotype IgG1, kappa. TC665 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 50 tests (20 µL of conjugate for up to 10<sup>6</sup> leucocytes from normal human peripheral blood).

**Specificity**

Anti-CD103, Ber-ACT8, was included in the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD103 (8).

Anti-CD11c, KB90, was included in the Sixth International Workshop and Conference on Leucocyte Cell Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD11c (9). Anti-CD11c, KB90, is reactive to the C-terminal region of CD11c (10).

Anti-CD19, HD37, was included in the Second, Third, Fourth and Fifth International Workshops and Conferences on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD19 (11, 12). Anti-CD19, HD37, labels human B cells in peripheral blood, bone marrow and other tissues. Additionally, B-cell lymphoproliferative disorders gave positive reactions with the antibody, i.e. acute lymphoblastic leukaemia, chronic lymphocytic leukaemia, hairy cell leukaemia, lymphoblastic lymphoma (Burkitt type), centroblastic/centrocytic lymphoma, and diffuse non-Hodgkin's lymphoma (11, 12).

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage**

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.

**Staining procedure**

1. Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
2. Add 20 µL of TC665 and mix gently by using a vortex mixer.
3. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.
4. Add 100 µL of Uti-Lyse™ Reagent A (Dako code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 1 mL of Uti-Lyse™ Reagent B (Dako code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
7. Add 2 mL of PBS (Dako code S3024) to the tube and resuspend the cells by using a vortex mixer.
8. Repeat step 6.
9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.
10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

**Procedural notes**

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the conjugated antibody reagent. The recommended Triple-Colour Control Reagent for TC665 is Dako code X0978.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, e.g. Dako EasyLyse™, code No. S2364, the PBS in step 9 should contain 1% paraformaldehyde or be analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis. The combination of FITC, RPE and APC is especially easy to work with since the emission spectra of these fluorochromes show minimum spectral overlap.

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Monoclonal Mouse Anti-Human CD19, clone HD37, conjugué à l'allophycocyanine (APC).

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
- van Dongen JJM, Adriaansen HJ. Immunobiology of leukemia. In: Henderson ES, Lister TA, Greaves MF, editors. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
- Kruschwitz M, Fritzsche G, Schwarting R, Micklem K, Mason DY, Falini B, et al. Ber-ACT8: New monoclonal antibody to the mucosa lymphocyte antigen. J Clin Pathol 1991;44:636-45.
- Sánchez ML, Almeida J, Vidrales B, López-Berges MC, García-Marcos MA, Moro MJ, et al. Incidence of phenotypic aberrations in a series of 467 patients with B chronic lymphoproliferative disorders: basis for the design of specific four-color stainings to be used for minimal residual disease investigation. Leukemia 2002;16:1460-9.
- Marotta G, Raspadori D, Sestigiani C, Scalia G, Bigazzi C, Lauria F. Expression of the CD11c antigen in B-cell chronic lymphoproliferative disorders. Leuk Lymphoma 2000;37:145-9.
- Pezzutto A, Dörken B, Rabinovitch PS, Ledbetter JA, Moldenhauer G, Clark EA. CD19 monoclonal antibody HD37 inhibits anti immunoglobulin-induced B cell activation and proliferation. J Immunol 1987;138:2793-9.
- Sato S, Tedder TF. BC3. CD19 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 133-5.
- Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematology neoplasias: results of an international consensus meeting. Cytometry 2001;46:23-7.
- Cepek KL, Wong DA, Brenner MB, Springer TA. AS7/8.4. CD103 (αE) cluster report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 1666-7.
- Hogg N. AS4. CD11c workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 347-8.
- Luk J, Luther E, Diamond MS, Springer TA. AS5.9. Subunit specificity and epitope mapping of Mac-1 and p150,95 mAb using chimeric CD11b X CD11c transfectants. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 2. p. 1599-1601.
- Nadler LM. B cell/leukemia panel workshop: summary and comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leukocyte Differentiation Antigens; 1984 Sept 17-20; Boston, USA. New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 3-43.
- Mason DY, Ladyman H, Gatter KC. Immunohistochemical analysis of monoclonal anti-B cell antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leukocyte Differentiation Antigens; 1984 Sept 17-20; Boston, USA. New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 245-55.

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Agilent Technologies Singapore (International) Pte Ltd.