MultiMix™
Triple-Colour Reagent
Anti-Human CD2/FITC
Anti-Human CD34/RPE
Anti-Human CD5/APC
Code No./ Ref./ Code-Nr. TC666

ENGLISH

Intended use

For in vitro diagnostic use.

TC666 is intended for use in flow cytometry. The reagent is intended for use in identification of cells expressing CD2, CD34 and CD5. CD2 is a pan-T-cell antigen and is present on virtually all normal T lymphocytes. CD2 is a useful marker in the assessment of lymphoid malignancies as it is expressed in the majority of precursor and postthymic lymphomas and leukaemias. In some neoplastic T-cell populations, e.g. in peripheral T-cell lymphomas, CD2 may be aberrantly deleted (1). CD34 is expressed by early haematopoietic stem and progenitor cells and also on endothelial cells. Antibodies to CD34 can be utilized to quantitate and purify haematopoietic stem/progenitor cells for transplantation and research (2-7). CD34 is a useful marker for stem and progenitor cells as well as for subclassification of leukaemias (3, 5, 6). CD5 is expressed by all mature T cells and almost all thymocytes. It is also present on a small subpopulation of normal and neoplastic B cells. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Reagent provided

TC666 comprises the following three matched fluorescent antibodies:

Monoclonal Mouse Anti-Human CD2, Clone MT910, conjugated with fluorescein isothiocyanate isomer 1 (FITC).
 Monoclonal Mouse Anti-Human CD34, Clone BIRMA-K3, conjugated with R-phycoerythrin (RPE).
 Monoclonal Mouse Anti-Human CD5, Clone DK23, conjugated with allophycocyanin (APC).

The three conjugates in TC666 have been produced from purified monoclonal mouse antibodies of isotype IgG1, kappa. TC666 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 50 tests (20 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).

Specificity

Anti-CD2, MT910, was included in the Second, Third and Fourth International Workshops and Conferences on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD2 (8).

Anti-CD34, BIRMA-K3, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by different laboratories confirmed its reactivity with CD34 and the class III epitope (2).

When tested by flow cytometry on buffy coat cells of peripheral blood as well as on gated lymphocyte populations, Anti-CD5, DK23, labels T cells. The specificity of Anti-CD5, DK23, is equivalent to that of the CD5-clustered antibody Leu-1 (9).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.

Staining procedure

1. Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
2. Add 20 µL of TC666 and mix gently by using a vortex mixer.
3. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.
4. Add 100 µL of Uti-Lyse™ Reagent A (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 1 mL of Uti-Lyse™ Reagent B (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
7. Add 2 mL of PBS (Dako code No. S3024) to the tube and resuspend the cells by using a vortex mixer.
8. Repeat step 6.
9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.
10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the isotypes and fluorochromes of the conjugated antibody reagent. The recommended Triple-Colour Control Reagent for TC666 is Dako code No. X0978.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as Dako EasyLyse™, code No. S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis. The combination of FITC, RPE and APC is especially easy to work with since the emission spectra of these fluorochromes show minimum spectral overlap.

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
1. Leong A, Cooper K, Leong F. Manual of diagnostic antibodies for immunohistology; London: Oxford University Press; 1999. p. 45-6.
2. Nishio H, Tada J, Hashiyama M, Hirn J, Ingles-Estevé J, Suda T. MC7. CD34 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 974-84.
3. Civin CI, Trischmann TM, Fackler MJ, Bernstein ID, Bühring HJ, Campos L, et al. M7.1. Report on the CD34 cluster workshop. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 818-25.
4. Höffkes H-G, Lowe JA, Pedersen RO, Schmidte G, McDonald DF. BIRMA-K3, a new monoclonal antibody for CD34 immunophenotyping and stem and progenitor cell assay. J Hematother 1996;5:261-70.
5. Civin CI, Strauss LC, Fackler MJ, Trischmann TM, Wiley JM, Loken MR. Positive stem cell selection - basic science. Prog Clin Biol Res 1990;333:387-402.
6. Krause DS, Fackler MJ, Civin CI, May WS. CD34: structure, biology, and clinical utility (review). Blood 1996;87:1-13.
7. Campos L, Guyotat D, Archimbaud E, Devaux Y, Treille D, Larese A, et al. Surface marker expression in adult acute myeloid leukaemia: correlations with initial characteristics, morphology and response to therapy. Br J Haematol 1989;72:161-6.
8. Dörken B, Möller P, Pezzutto A, Schwartz-Albiez R, Moldenhauer G. B3. B-cell antigens: CD20. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p 46-8.
9. Stein H, Lennert K, Feller AC, Mason DY. Immunohistological analysis of human lymphoma: correlation of histological and immunological categories. Adv Cancer Res 1984;42:67-147.

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Agilent Technologies Singapore (International) Pte Ltd.