**MultiMix™****Triple-Colour Reagent****Anti-Human Plasma Cell/FITC****Anti-Human Lambda Light Chains/RPE****Anti-Human Kappa Light Chains/APC****Code No./ Ref./ Code-Nr. TC670****ENGLISH****Intended use**

For in vitro diagnostic use.

TC670 is intended for use in flow cytometry. The reagent is intended for intracellular identification of the plasma cell antigen, p63 protein, and kappa and lambda light chain expression. Anti-Human Plasma Cell, VS38c, labels normal and neoplastic plasma cells (1, 2). The kappa and lambda light chain antibodies allow evaluation of intracellular light chain expression in VS38c-positive plasma cells, which is of value in detection of clonality in neoplastic plasma cells. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Reagent provided**

TC670 comprises the following three matched fluorescent antibodies:

Monoclonal Mouse Anti-Human Plasma Cell, Clone VS38c, conjugated with fluorescein isothiocyanate isomer 1 (FITC).

Polyclonal Rabbit Anti-Human Lambda Light Chains, conjugated with R-phycoerythrin (RPE).

Polyclonal Rabbit Anti-Human Kappa Light Chains, conjugated with allophycocyanin (APC).

The FITC conjugate in TC670 has been produced from a purified monoclonal mouse antibody of isotype IgG1, kappa. The RPE and APC conjugates in TC670, Anti-Lambda Light Chains and Anti-Kappa Light Chains, have been produced from F(ab')<sub>2</sub> fragments of affinity-isolated polyclonal rabbit antibodies. TC670 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaNa<sub>3</sub>, pH 7.2. Each vial contains 50 tests (20 µL of conjugate for up to 10<sup>6</sup> cells (a mixture of the cell lines HS-Sultan, U266 and Amo-1)).

**Immunogens**

Monoclonal Mouse Anti-Human Plasma Cell, Clone VS38c: Membrane preparation of MCF-7 AdR cells (breast carcinoma cell line) (2).

Rabbit Anti-Human Kappa Light Chains: Polyclonal immunoglobulin light chains of kappa type isolated from a pool of human sera.

Rabbit Anti-Human Lambda Light Chains: Polyclonal immunoglobulin light chains of lambda type isolated from a pool of human sera.

**Specificity**

Anti-Plasma Cell, VS38c, was submitted to and characterized as anti-p63 protein at the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens (1).

Anti-Lambda Light Chains reacts with free lambda chains as well as lambda chains in intact immunoglobulin molecules. The specificity has been ascertained as described below. To obtain maximum sensitivity, the crossed immunoelectrophoresis specificity test was performed prior to affinity purification and F(ab')<sub>2</sub> fragmentation.

Crossed immunoelectrophoresis: Only lambda-related precipitates appear when the antibody is tested against human plasma. Staining: Coomassie Brilliant Blue.

Flow cytometry: When Anti-Lambda Light Chains is applied as described in the enclosed staining procedure in combination with Anti-CD19/RPE, HD37, on lysed human whole blood, a specific staining of a part of the CD19-positive B lymphocytes is seen corresponding to the expected range of lambda light chain expression.

Anti-Kappa Light Chains reacts with free kappa chains as well as kappa chains in intact immunoglobulin molecules. The specificity has been ascertained as described below. To obtain maximum sensitivity, the crossed immunoelectrophoresis specificity test was performed prior to affinity purification and F(ab')<sub>2</sub> fragmentation.

Crossed immunoelectrophoresis: Only kappa-related precipitates appear when the antibody is tested against human plasma. Staining: Coomassie Brilliant Blue.

Flow cytometry: When Anti-Kappa Light Chains is applied as described in the enclosed staining procedure in combination with Anti-CD19/RPE, HD37, on lysed human whole blood, a specific staining of a part of the CD19-positive B lymphocytes is seen corresponding to the expected range of kappa light chain expression.

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage**

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

Before staining samples of peripheral blood, the mononuclear cells must be isolated by centrifugation on a separation medium or the blood sample washed to remove the soluble serum proteins. As human monocytes bind serum immunoglobulins via their surface Fc receptors, these cells should be removed by depletion or identified.

1. Transfer 50 µL (up to 10<sup>6</sup> cells) anticoagulated (EDTA) blood cells to a 12 mm x 75 mm polystyrene test tube.
2. Add 2 mL of PBS (Dako code No. S3024). Mix gently by using a vortex mixer.
3. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
4. Repeat steps 2 and 3 two more times.
5. Add 10 µL Rabbit Immunoglobulin Fraction (Dako code No. X0903) for blocking. Mix gently by using a vortex mixer and incubate the tube at 37 °C for 30 minutes.
6. Add 100 µL IntraStain Reagent A (Fixation) (Dako code No. K2311) to the tube and mix gently by using a vortex mixer to ensure that the cells are in suspension. Incubate the tube for 15 minutes at room temperature.
7. Add 2 mL PBS (Dako code No. S3024) and mix gently by using a vortex mixer.
8. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
9. Mix thoroughly by using a vortex mixer to ensure that the cells are in suspension.
10. Add 100 µL IntraStain Reagent B (Permeabilization) (Dako code No. K2311) to the tube.
11. Then add 20 µL of TC670 and mix gently by using a vortex mixer.
12. Incubate the tube for 15 minutes at room temperature in the dark.
13. Repeat steps 7 and 8.
14. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.

15. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 8 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

#### Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 11: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the conjugated antibody reagent. The recommended Triple-Colour Control Reagent for TC670 is Dako code No. X0979.

Step 15: Abnormal numbers of cells expressing the target antigen or cell showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis. The combination of FITC, RPE and APC is especially easy to work with since the emission spectra of these fluorochromes show minimum spectral overlap.

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Agilent Technologies Singapore (International) Pte Ltd.