MultiMix™
Triple-Colour Reagent
Anti-Human CD38/FITC
Anti-Human CD56/RPE
Anti-Human CD45/APC
Code TC671

ENGLISH

Intended use	For in vitro diagnostic use. TC671 is intended for use in flow cytometry. The reagent allows simultaneous detection and enumeration of CD38-positive and CD56-positive cells in combination with CD45-positive leucocytes. CD38 is a useful marker in the immunophenotyping of acute leukaemias. Additionally, antibodies to CD38 are valuable for the identification of plasma cells, as poorly differentiated plasma cells may mimick other blastic lymphoid cells (1). CD56 is the prototypic marker of natural killer (NK) cells and it is also present on a subset of CD4-positive and CD8-positive T cells in peripheral blood (2). CD45 is one of the most abundant leucocyte cell surface glycoproteins and is expressed exclusively on cells of the haematopoietic system and their progenitors (3). Anti-CD45, together with a panel of other antibodies, is considered essential for the initial evaluation of chronic lymphoproliferative disorders and acute leukaemias (4). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Reagent provided	TC671 comprises the following three matched fluorescent antibodies: Monoclonal Mouse Anti-Human CD38, Clone AT13/5, conjugated with fluorescein isothiocyanate isomer 1 (FITC). Monoclonal Mouse Anti-Human CD56, Clone C5.9, conjugated with R-phycoerythrin (RPE). Monoclonal Mouse Anti-Human CD45, Clone 2D1, conjugated with allophycocyanin (APC). The three conjugates in TC671 have been produced from purified monoclonal mouse antibodies of isotypes IgG1, kappa (Anti-CD38), IgG2b, kappa (Anti-CD56) and IgG1, kappa (Anti-CD45). TC671 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 50 tests (20 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).
Specificity	Anti-CD38, AT13/5, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD38 (5). Anti-CD56, C5.9, labels NK cells in flow cytometry when analysing buffy coat cells of peripheral blood as well as gated lymphocyte populations and NK cell cultures. The staining pattern parallels that observed with Anti-CD56, MOC-1, which has been clustered as an antibody to CD56 (6). Anti-CD45, 2D1, was included in the Third International Workshop and Conference on Human Leucocyte Differentiation Antigens (7).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact Dako Technical Support.
Staining procedure	1. Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube. 2. Add 20 µL of TC671 and mix gently by using a vortex mixer. 3. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes. 4. Add 100 µL of Uti-Lyse™ Reagent A (Dako Code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark. 5. Add 1 mL of Uti-Lyse™ Reagent B (Dako Code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark. 6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube. 7. Add 2 mL of PBS (Dako Code S3024) to the tube and resuspend the cells by using a vortex mixer. 8. Repeat step 6. 9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS. 10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis. Procedural notes Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control. Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory. It is optional to include a control reagent test tube. The control reagent should match the conjugated antibody reagent.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, e.g. Dako EasyLyse™, Code S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis. The combination of FITC, RPE and APC is especially easy to work with since the emission spectra of these fluorochromes show minimum spectral overlap.

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1. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 87-8.
2. Poggi A. CD56. CD Guide. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 805-6.
3. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 95-8.
4. Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematolymphoid neoplasias: results of an international consensus meeting. Cytometry 2001;46:23-7.
5. Morra M, Deaglio S, Mallone R, di Rosa G, Tibaldi E, Zini G, et al. BC10. CD38 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 151-4.
6. Costa P. NK4. CD56 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 271-2.
7. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 788-803.

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Agilent Technologies Singapore (International) Pte Ltd.