MultiMix™**Triple-Colour Reagent****Anti-Human B Cell (FMC7)/FITC****Anti-Human CD23/RPE****Anti-Human CD19/APC****Code TC683****ENGLISH****Intended use**

For in vitro diagnostic use.

TC683 is intended for use in flow cytometry. The reagent is intended for use in identification of FMC7-, CD23- and CD19-positive cells.

The FMC7 antigen is present on a subpopulation of relatively mature B cells in normal peripheral blood (1, 2). The antigen is expressed in a range of lymphoproliferative disorders, including B-cell prolymphocytic leukemia (B-PLL), hairy cell leukemia (HCL), mantle cell lymphoma (MCL), follicle centre lymphoma (FCL), marginal zone lymphoma (MZL), mucosa-associated lymphoid tissue (MALT) lymphoma, diffuse large B cell lymphoma (DLCL), Burkitt's lymphoma (BL) and Waldenström macroglobulinaemia (WM) (1-4). FMC7 has a low expression or is absent in B-cell chronic lymphocytic/small lymphocytic lymphoma (B-CLL/SL), null-cell acute lymphoblastic leukemia (Null-ALL), B-cell acute lymphoblastic leukemia (B-ALL) and multiple myeloma (MM) (1-5). The FMC7 antigen is not expressed in acute myeloid leukemia (AML), chronic myeloid leukemia (CML), T- ALL and T-PLL (5).

CD23 is identical to the low affinity IgE receptor (FcεRII) found on B cells. CD23 is primarily expressed on B cells and monocytes and shows a strong expression on EBV-transformed B lymphoblasts (6). CD23 is also present on a large variety of other cells, such as T cells, eosinophils, platelets, Langerhans cells and a subset of thymic epithelial cells (6). CD23 is expressed on neoplastic cells from cases of B-cell chronic lymphocytic leukemia and some cases of centroblastic/centrocytic lymphoma but is not found in other types of lymphoid neoplasms (7).

CD19 is the broadest lineage-specific surface marker for B cells and it is present on the surface of virtually all B lymphocytes, including early B progenitor cells. CD19 expression is maintained in B-lineage cells that have undergone neoplastic transformation (8). Antibodies to CD19 are considered essential for the initial evaluation of acute and chronic lymphoproliferative disorders (9).

Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Reagent provided

TC683 comprises the following three matched fluorescent antibodies:

Monoclonal Mouse Anti-Human B Cell (FMC7), Clone FMC7, conjugated with fluorescein isothiocyanate isomer 1 (FITC).

Monoclonal Mouse Anti-Human CD23, Clone MHM6, conjugated with R-phycoerythrin (RPE).

Monoclonal Mouse Anti-Human CD19, Clone HD37, conjugated with allophycocyanin (APC).

The three conjugates in TC683 have been produced from purified monoclonal mouse antibodies of isotypes IgM, kappa (Anti-B Cell (FMC7)), IgG1, kappa (Anti-CD23) and IgG1, kappa (Anti-CD19). TC683 is provided in liquid form in buffer containing 0.4% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 8.0. Each vial contains 50 tests (20 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).

Specificity

Anti-B Cell (FMC7) reacts with 3-6% of normal peripheral blood mononuclear cells. These represent a subset of B cells that also expresses surface membrane immunoglobulin (Smlg). The antibody does not label B cells without Smlg, T cells and monocytes (5). In lymphoproliferative disorders, Anti-B Cell (FMC7) has been shown to label B-PLL (2/2 cases (1)), HCL (14/14 cases (1)), MCL (15/15 cases (1) and 12/13 cases (2)), FCL (24/24 cases (1) and 46/62 cases (2)), MZL (20/20 cases (1) and 13/14 cases (2)), MALT lymphoma (14/17 cases (2)), DLCL (54/64 cases (2)), BL (2/2 cases (1) and 1/1 case (2)) and WM (5/6 cases (1) and 10/11 cases (4)).

Anti-B Cell (FMC7) shows weak or no reactivity in B-CLL/SSL (27/27 cases (2), 106/121 cases (1); and 17/17 cases (5)), Null-ALL (19/19 cases (5)), B-ALL (1/1 case (2)), T-ALL (5/5 cases (5)), AML/CML (6/6 cases (5)), T-PLL (1/1 case (5)) and MM (26/26 cases (4)).

Anti-CD23, MHM6, was included in the Second, Third, Fourth and Fifth International Workshops and Conferences on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD23 (6).

Anti-CD19, HD37, was included in the Second, Third, Fourth and Fifth International Workshops and Conferences on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD19 (10,11).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact Dako Technical Support.

Staining procedure

1. Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
2. Add 20 µL of TC683 and mix gently by using a vortex mixer.
3. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.
4. Add 100 µL of Uti-Lyse™ Reagent A (Dako Code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 1 mL of Uti-Lyse™ Reagent B (Dako Code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
7. Add 2 mL of PBS (Dako Code S3024) to the tube and resuspend the cells by using a vortex mixer.
8. Repeat step 6.
9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.
10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the isotypes and fluorochromes of the conjugated antibody reagent.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as Dako EasyLyse™, Code S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolor reagents are preferable to single-color reagents for the comprehensive analysis of flow cytometry specimens. The correct use of color compensation is particularly important in multicolor analysis. The combination of FITC, RPE and APC is especially easy to work with since the emission spectra of these fluorochromes show minimum spectral overlap.

Product-specific limitations

Lysing of blood samples with lysing reagents containing fixative prior to incubation with TC683 can destroy the epitope recognized by Anti-Human B Cell (FMC7)/FITC.

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1. Ahmad E, Garcia D, Davis BH. Clinical utility of CD23 and FMC7 antigen coexistent expression in B-cell lymphoproliferative disorder subclassification. *Cytometry* 2002;50:1-7.
2. Garcia DP, Rooney MT, Ahmad E, Davis BH. Diagnostic usefulness of CD23 and FMC-7 antigen expression patterns in B-cell lymphoma classification. *Am J Clin Pathol* 2001;115:258-65.
3. Zola H, Neoh SH, Potter A, Melo JV, De Oliveria MSP, Catovsky D. Markers of differentiated B cell leukaemia: CD22 antibodies and FMC7 react with different molecules. *Dis Markers* 1987;5:227-35.
4. San Miguel JF, Caballero MD, Gonzalez M, Zola H, Lopez Borrascas A. Immunological phenotype of neoplasms involving the B cell in the last step of differentiation. *Br J Haematol* 1986;62:75-83.
5. Brooks DA, Beckman IGR, Bradley J, McNamara PJ, Thomas ME, Zola H. Human lymphocyte markers defined by antibodies derived from somatic cell hybrids. IV. A monoclonal antibody reacting specifically with a subpopulation of human B lymphocytes. *J Immunol* 1981;126:1373-7.
6. Sarfati M, Ishihara H, Delespesse G. B7. CD23 workshop panel report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leucocyte typing V. White celldifferentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA.* Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 1. p. 530-3.
7. Pallesen G. B2.5. The distribution of CD23 in normal human tissues and in malignant lymphomas. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. *Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England.* Oxford, New York, Tokyo: Oxford University Press; 1987. p. 383-6.
8. Sato S, Tedder TF. BC3. CD19 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan.* New York, London: Garland Publishing Inc.; 1997. p. 133-5.
9. Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematology neoplasias: results of an international consensus meeting. *Cytometry* 2001;46:23-7.
10. Nadler LM. B cell/leukemia panel workshop: summary and comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. *Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leukocyte Differentiation Antigens; 1984 Sept 17-20; Boston, USA.* New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 3-43.
11. Mason DY, Ladyman H, Gatter KC. Immunohistochemical analysis of monoclonal anti-B cell antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. *Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leukocyte Differentiation Antigens; 1984 Sept 17-20; Boston, USA.* New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 245-55.

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Agilent Technologies Singapore (International) Pte Ltd.