MultiMix™

Triple-Colour Reagent
Anti-Human CD13/FITC
Anti-Human HLA-DR Antigen/RPE
Anti-Human CD117/APC
Code TC685

ENGLISH**Intended use**

For in vitro diagnostic use.

TC685 is intended for use in flow cytometry. The reagent is intended for use in identification of cells expressing CD13, HLA-DR Antigen and CD117.

CD13 is expressed on the surface of committed granulocyte-monocyte progenitors (CFU-GM) and by cells of the granulocyte and monocyte lineages at all stages of differentiation, as well as by the neoplastic counterparts of these cells (1, 2). The antigen is also expressed on a small proportion of large granular lymphocytes but not on other lymphocytes. Further, CD13 is found in some non-haematopoietic tissues (1-4).

The HLA-DR antigen is constitutively expressed on antigen-presenting cells, such as B lymphocytes, monocytes and dendritic cells but can also be detected on activated T lymphocytes and activated granulocytes (5, 6). Occasionally, natural killer cells express HLA-DR antigen (1). The antigen has been found expressed in cases of different types of acute lymphoblastic leukaemias, acute myeloid leukaemias (AML) except AML-M3, chronic T-cell leukaemias, chronic myeloid leukaemias (CML) and B- and T-cell non-Hodgkin's leukaemias (1, 7).

CD117 is useful for identification of AML, and for classification of leukaemias (8-10). CD117 is a marker for tissue mast cells, haematopoietic stem cells, and progenitor cells in normal human bone marrow. The majority (50-70%) of CD117-positive marrow cells coexpress CD34 and comprise progenitor cells and their precursors of all haematopoietic lineages (10).

TC685 is not intended for tissue typing. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Reagent provided

TC685 comprises the following three matched fluorescent antibodies:

Monoclonal Mouse Anti-Human CD13, Clone WM-47, conjugated with fluorescein isothiocyanate isomer 1 (FITC).

Monoclonal Mouse Anti-Human HLA-DR Antigen, Clone AB3, conjugated with R-phycoerythrin (RPE).

Monoclonal Mouse Anti-Human CD117, Clone 104D2, conjugated with allophycocyanin (APC).

The three conjugates in TC685 have been produced from purified monoclonal mouse antibodies of isotypes IgG1, kappa (Anti-CD13), IgG2a, kappa (Anti-HLA-DR Antigen) and IgG1, kappa (Anti-CD117). TC685 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 50 tests (20 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).

Specificity

Anti-CD13, WM-47, was included in the Fourth and Fifth International Workshops and Conferences on Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with the CD13 antigen (4, 11). Anti-CD13, WM-47, reacts with CFU-GM cells, normal granulocytic and monocytic cells at all stages of differentiation (1, 4, 11). It does not react with lymphocytes and platelets (3). Flow cytometric analysis demonstrated that Anti-CD13, WM-47, labels AML (9/9 cases) and CML in myeloid blast crisis (7/7 cases). It does not react with common acute lymphoblastic leukaemia (0/8 cases) or T-cell acute lymphoblastic leukaemia (0/1 case) (3).

Anti-Human HLA-DR Antigen, AB3, reacts with the DRw52 determinant of the HLA-DR molecule. The specificity has been tested with transfectants expressing HLA-DRw52. Anti-Human HLA-DR Antigen, AB3, does not react with transfectants expressing HLA-DP or HLA-DQ antigen. The antibody has been submitted to the Eighth International Workshop and Conference on Human Leucocyte Differentiation Antigens.

Anti-CD117, 104D2, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD117 (10).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.

Staining procedure

1. Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
2. Add 20 µL of TC685 and mix gently by using a vortex mixer.
3. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.
4. Add 100 µL of Uti-Lyse™ Reagent A (Dako code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 1 mL of Uti-Lyse™ Reagent B (Dako code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
7. Add 2 mL of PBS (Dako code S3024) to the tube and resuspend the cells by using a vortex mixer.
8. Repeat step 6.
9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.
10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the isotypes and fluorochromes of the conjugated antibody reagent. The recommended Triple-Colour Control Reagent for TC685 is Dako code X7904.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as Dako EasyLyse™, code S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis. The combination of FITC, RPE and APC is especially easy to work with since the emission spectra of these fluorochromes show minimum spectral overlap.

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
1. van Dongen JJM, Adriaansen HJ. Immunobiology of leukemia. In: Henderson ES, Lister TA, Greaves MF, editors. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
2. Turni L, Shaw S, Watson B, Mason D. CD guide. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 June 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 761.
3. Favaloro EJ, Bradstock KF, Kabral A, Grimsley P, Zowtyj H, Zola H. Further characterization of human myeloid antigens (gp160,95; gp150; gp67): investigation of epitopic heterogeneity and non-haemopoietic distribution using panels of monoclonal antibodies belonging to CD-11b, CD-13 and CD-33. *Br J Haematol* 1988;69:163-71.
4. Look AT, Ashmun RA, Shapiro LH, O'Connell PJ, Gerkis V, D'Apice AJ, et al. M2.1. Report on the CD13 (aminopeptidase N) cluster workshop. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 784-7.
5. Naverrete CV. The HLA system in blood transfusion. *Baillière's Clin Haematol* 2000;13:511-32.
6. Erber WN, Pinching AJ, Mason DY. Immunocytochemical detection of T and B cell populations in routine blood smears. *Lancet* 1984;i:1042-5.
7. Smith MEF, Holgate CS, Williamson JMS, Grigor I, Quirke P, Bird CC. Major histocompatibility complex class II antigen expression in B and T cell non-Hodgkin's lymphoma. *J Clin Pathol* 1987;40:34-41.
8. Bahia DMM, Yamamoto M, de Lourdes M, Chauffaille LF, Kimura EYS, Bordin JO, et al. Aberrant phenotypes in acute myeloid leukaemia: a high frequency and clinical significance. *Haematologica* 2001;86:801-6.
9. Bühring H-J, Ullrich A, Schaudt K, Müller CA, Busch FW. The product of the proto-oncogene c-kit (P145^{c-kit}) is a human bone marrow surface antigen of hemopoietic precursor cells which is expressed on a subset of acute nonlymphoblastic leukemic cells. *Leukemia* 1991;5:854-60.
10. Ashman LK, Cambareri AC, Nguyen L, Bühring H-J. CR5. CD117 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 816-8.
11. Ashmun RA, Holmes KV, Shapiro LH, Favaloro EJ, Razak K, de Crom RPG, et al. M3. CD13 (aminopeptidase N) cluster workshop report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 1. p. 771-5.

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Agilent Technologies Singapore (International) Pte Ltd.