MultiMix™**Triple-Colour Reagent****Anti-Human CD41/FITC****Anti-Human CD34/RPE****Anti-Human CD61/APC****Code No./ Code/ Code-Nr. TC687****ENGLISH****Intended use**

For in vitro diagnostic use.

TC687 is intended for use in flow cytometry. The reagent is intended for identification of CD41, CD61 and CD34 positive cells.

CD41 and CD61 are selective markers of platelets and platelet precursors, and they may be of value for immunophenotyping of megakaryoblastic leukaemias (1, 2). The CD41/CD61 complex appears early in megakaryocyte maturation (3). The activated CD41/CD61 complex is a receptor for von Willebrand factor, soluble fibrinogen and fibronectin and plays a central role in platelet activation and aggregation (1, 2, 4).

The CD34 antigen is present on human haematopoietic progenitor cells in bone marrow and blood. In contrast, the CD34 antigen is normally not detected on peripheral blood leucocytes and platelets (5). Approximately 60% of acute B-lymphoid leukaemias and 40% of acute myeloid leukaemias (AML), and 1% to 5% of acute T-lymphoid leukaemias express CD34 (5). Chronic lymphoid leukaemias, lymphomas and multiple myelomas have been found to be uniformly CD34 negative (5, 6).

Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Reagent provided

TC687 comprises the following three matched fluorescent antibodies:

Monoclonal Mouse Anti-Human CD41, Clone 5B12, conjugated with fluorescein isothiocyanate isomer 1 (FITC).

Monoclonal Mouse Anti-Human CD34, Clone BIRMA-K3, conjugated with R-phycoerythrin (RPE).

Monoclonal Mouse Anti-Human CD61, Platelet Glycoprotein IIIa, Clone Y2/51, conjugated with allophycocyanin (APC).

The three conjugates in TC687 have been produced from purified monoclonal mouse antibodies of isotypes IgG1, kappa (Anti-CD41), IgG1, kappa (Anti-CD34) and IgG1, kappa (Anti-CD61). TC687 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 50 tests (20 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).

Specificity

Anti-CD41, 5B12, was included in the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD41. Erroneously the clone was listed as SB12 (7).

Anti-CD34, BIRMA-K3, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by different laboratories confirmed its reactivity with CD34 and the class III epitope (8). In flow cytometry the antibody labels KG-1a cells (a primitive myeloid leukaemia cell line, known to be CD34-positive) (9).

Anti-CD61, Y2/51, was included in the Fourth International Workshop and Conference on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD61 (10).

Anti-CD61 detects platelets in peripheral blood and bone marrow and it also reacts with megakaryocytes (11).

Precautions

1. For professional users.

2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.

3. As with any product derived from biological sources, proper handling procedures should be used.

Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.

Staining procedure

1. Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
2. Add 20 µL of TC687 and mix gently by using a vortex mixer.
3. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.
4. Add 100 µL of Uti-Lyse™ Reagent A (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 1 mL of Uti-Lyse™ Reagent B (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
7. Add 2 mL of PBS (Dako code No. S3024) to the tube and resuspend the cells by using a vortex mixer.
8. Repeat step 6.
9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.
10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the isotypes and fluorochromes of the conjugated antibody reagent. The recommended Triple-Colour Control Reagent for TC687 is Dako code No. X0978.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as Dako EasyLyse™, code No. S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis. The combination of FITC, RPE and APC is especially easy to work with since the emission spectra of these fluorochromes show minimum spectral overlap.

Product-specific limitations

TC687 will in conjunction with Dako lysing procedures form platelet-leukocyte aggregates (PLA). High levels of PLA may be seen with Dako EasyLyse™ lyse-wash procedure. PLA can be reduced by lysing prior to addition of TC687 or by immediately fixing the sample and then employing a no-lyse-no-wash procedure (12).

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1. van Dongen JJM, Adriaansens HJ. Immunobiology of leukemia. In: Henderson ES, Lister TA, Greaves MF, editors. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
2. Sun QH, Newman PJ. CD Guide. CD61. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 810.
3. Vinci G, Tabilio A, Deschamps JF, van Haeke D, Henri A, Guichard J, et al. Immunological study of *in vitro* maturation of human megakaryocytes. Br J Haematol 1984;56:589-605.
4. Nachman RL, Leung LLK. Complex formation of platelet membrane glycoproteins IIb and IIIa with fibrinogen. J Clin Invest 1982;69:263-9.
5. Civin CI, Trischmann TM, Fackler MJ, Bernstein ID, Bühring HJ, Campos L, et al. M7.1. Report on the CD34 cluster workshop. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 818-25.
6. Campos L, Guyotat D, Archimbaud E, Devaux Y, Treille D, Larese A, et al. Surface marker expression in adult acute myeloid leukaemia: correlations with initial characteristics, morphology and response to therapy. Br J Haematol 1989;72:161-6.
7. Honda S, Felding-Habermann B, Loftus J, Annis D, Kunicki TJ. CD41/CD61 cluster workshop report: localization of epitopes on integrins $\alpha_{IIb}\beta_3$ (CD41/CD61) and $\alpha_v\beta_3$ (CD51/CD61). In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 1293-8.
8. Nishio H, Tada J, Hashiyama M, Hirn J, Ingles-Esteve J, Suda T. MC7. CD34 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 974-84.
9. Höffkes H-G, Lowe JA, Pedersen RO, Schmidtko G, McDonald DF. BIRMA-K3, a new monoclonal antibody for CD34 immunophenotyping and stem and progenitor cell assay. J Hematother 1996;5:261-70.
10. von dem Borne AEG, Modderman PW, Admiraal LG, Nieuwenhuis HK. Joint report of the platelet section. P1: Platelets antibodies, the overall results. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 951-66.
11. Gatter KC, Cordell JL, Turley H, Heryet A, Kieffer N, Anstee DJ, et al. The immunohistological detection of platelets, megakaryocytes and thrombi in routinely processed specimens. Histopathology 1988;13:257-67.
12. Hagberg IA, Lyberg T. Evaluation of circulating platelet-leukocyte conjugates: a sensitive flow cytometric assay well suited for clinical studies. Platelets 2000;11:151-60.

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Agilent Technologies Singapore (International) Pte Ltd.