MultiMix™**Triple-Colour Reagent
Anti-Human CD19/FITC
Anti-Human CD34/RPE
Anti-Human CD22/APC
Code No./ Code/ Code-Nr. TC689****ENGLISH****Intended use**

For in vitro diagnostic use.

TC689 is intended for use in flow cytometry. The reagent is intended for use in identification of CD19- and CD22-positive B cells and CD34-positive progenitor cells.

CD19 is the broadest lineage-specific surface marker for B cells and it is present on the surface of virtually all B lymphocytes, including early B progenitor cells (1). CD19 expression is maintained in B-lineage cells that have undergone neoplastic transformation (2). Antibodies to CD19 are considered essential for the initial evaluation of acute and chronic lymphoproliferative disorders (3).

The CD34 antigen is present on human haematopoietic progenitor cells in bone marrow and blood. In contrast, the CD34 antigen is normally not detected on peripheral blood leucocytes and platelets (4). Approximately 60% of acute B-lymphoid leukaemias and 40% of acute myeloid leukaemias (AML), and 1% to 5% of acute T-lymphoid leukaemias express CD34 (4). Chronic lymphoid leukaemias, lymphomas and multiple myelomas have been found to be uniformly CD34 negative (4, 5).

The CD22 antigen is expressed on normal and neoplastic B-lymphocytes in samples of human peripheral blood and bone marrow. CD22 is present in neoplasms of B cell origin, including primitive lymphomas, which lack monoclonal surface immunoglobulin, and in most cases of common acute lymphoblastic leukaemia (6).

Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Reagent provided

TC689 comprises the following three matched fluorescent antibodies:

Monoclonal Mouse Anti-Human CD19, Clone HD37, conjugated with fluorescein isothiocyanate isomer 1 (FITC).

Monoclonal Mouse Anti-Human CD34, Clone BIRMA-K3, conjugated with R-phycoerythrin (RPE).

Monoclonal Mouse Anti-Human CD22, Clone 4KB128, conjugated with allophycocyanin (APC).

The three conjugates in TC689 have been produced from purified monoclonal mouse antibodies of isotype IgG1, kappa. TC689 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L Na₃N, pH 7.2. Each vial contains 50 tests (20 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).

Specificity

Anti-CD19, HD37, was included in the Second, Third, Fourth and Fifth International Workshops and Conferences on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD19 (7, 8). Anti-CD19, HD37, labels human B cells in peripheral blood, bone marrow and other tissues. Additionally, B-cell lymphoproliferative disorders gave positive reactions with the antibody, i.e. acute lymphoblastic leukaemia, chronic lymphocytic leukaemia, hairy cell leukaemia, lymphoblastic lymphoma (Burkitt type), centroblastic/centrocytic lymphoma, and diffuse non-Hodgkin's lymphoma (7, 8).

Anti-CD34, BIRMA-K3, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD34 and the class III epitope (9). In flow cytometry the antibody labels KG-1a cells (a primitive myeloid leukaemia cell line, known to be CD34-positive) (10).

Anti-CD22, 4KB128, was included in the Third and Fifth International Workshops and Conferences on Human Leucocyte Differentiation Antigens (11, 12), and studies by a number of laboratories confirmed its reactivity with the CD22 antigen. Anti-CD22, 4KB128, reacts with normal and neoplastic B cells and does not react with other haematopoietic cell types (11).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.

Staining procedure

1. Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
2. Add 20 µL of TC689 and mix gently by using a vortex mixer.
3. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.
4. Add 100 µL of Uti-Lyse™ Reagent A (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 1 mL of Uti-Lyse™ Reagent B (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
7. Add 2 mL of PBS (Dako code No. S3024) to the tube and resuspend the cells by using a vortex mixer.
8. Repeat step 6.
9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.
10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the conjugated antibody reagent. The recommended Triple-Colour Control Reagent for TC689 is Dako code No. X0978.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, e.g. Dako EasyLyse™, code No. S2364, the PBS in step 9 should contain 1% paraformaldehyde or be analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis. The combination of FITC, RPE and APC is especially easy to work with since the emission spectra of these fluorochromes show minimum spectral overlap.

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
1. Pezzutto A, Dörken B, Rabinovitch PS, Ledbetter JA, Moldenhauer G, Clark EA. CD19 monoclonal antibody HD37 inhibits anti immunoglobulin-induced B cell activation and proliferation. *J Immunol* 1987;138:2793-9.
2. Sato S, Tedder TF. BC3. CD19 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference*; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 133-5.
3. Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematology neoplasias: results of an international consensus meeting. *Cytometry* 2001;46:23-7.
4. Civin CI, Trischmann TM, Fackler MJ, Bernstein ID, Bühring HJ, Campos L, et al. M7.1. Report on the CD34 cluster workshop. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. *Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference*; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 818-25.
5. Campos L, Guyotat D, Archimbaud E, Devaux Y, Treille D, Larese A, et al. Surface marker expression in adult acute myeloid leukaemia: correlations with initial characteristics, morphology and response to therapy. *Br J Haematol* 1989;72:161-6.
6. Mason DY, Stein H, Gerdes J, Pulford KAF, Ralfkiaer E, Falini B, et al. Value of monoclonal anti-CD22 (p135) antibodies for the detection of normal and neoplastic B lymphoid cells. *Blood* 1987;69:836-40.
7. Nadler LM. B cell/leukemia panel workshop: summary and comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. *Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leucocyte Differentiation Antigens*; 1984 Sept 17-20; Boston, USA. New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 3-43.
8. Mason DY, Ladyman H, Gatter KC. Immunohistochemical analysis of monoclonal anti-B cell antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. *Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leucocyte Differentiation Antigens*; 1984 Sept 17-20; Boston, USA. New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 245-55.
9. Nishio H, Tada J, Hashiyama M, Hirn J, Ingles-Esteve J, Suda T. MC7. CD34 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference*; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 974-84.
10. Höffkes H-G, Lowe JA, Pedersen RO, Schmidtko G, McDonald DF. BIRMA-K3, a new monoclonal antibody for CD34 immunophenotyping and stem and progenitor cell assay. *J Hematother* 1996;5:261-70.
11. Kehl JH. B6 CD22 Workshop Panel report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference*; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press. Volume 1. 1995. p. 523-5.
12. Moldenhauer G, Schwartz R, Dörken B, Hämmerling GJ. B2.2 Biochemical characterization and epitope analysis of B-lymphocyte specific surface antigens defined by clustering Workshop monoclonal antibodies. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. *Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference*; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 378-82.

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Agilent Technologies Singapore (International) Pte Ltd.