MultiMix™

Triple-Colour Reagent

Anti-Human CD3/FITC

Anti-Human CD19/RPE

Anti-Human CD45/APC

Code No./ Code/ Code-Nr. TC690

ENGLISH

Intended use

For in vitro diagnostic use.

TC690 is intended for use in flow cytometry. The reagent is intended for use in identification of CD3-positive T-cell subsets and CD19-positive B cells, in combination with CD45-positive leucocytes.

CD3 is a pan-T-cell restricted antigen, which is a valuable marker for normal and neoplastic T cells. Antibodies to CD3 are considered essential for the initial evaluation of acute and chronic lymphoproliferative disorders (1).

CD19 is the broadest lineage-specific surface marker for B cells and it is present on the surface of virtually all B lymphocytes, including early B progenitor cells (2). CD19 expression is maintained in B-lineage cells that have undergone neoplastic transformation (3). Antibodies to CD19 are considered essential for the initial evaluation of acute and chronic lymphoproliferative disorders (1).

CD45 is one of the most abundant leucocyte cell surface glycoproteins and is expressed exclusively on cells of the haematopoietic system and their progenitors (4). Anti-CD45, together with a panel of other antibodies, is considered essential for the initial evaluation of chronic lymphoproliferative disorders and acute leukaemias (1).

Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Reagent provided

TC690 comprises the following three matched fluorescent antibodies:

Monoclonal Mouse Anti-Human CD3, Clone UCHT1, conjugated with fluorescein isothiocyanate isomer 1 (FITC).

Monoclonal Mouse Anti-Human CD19, Clone HD37, conjugated with R-phycoerythrin (RPE).

Monoclonal Mouse Anti-Human CD45, Clone 2D1, conjugated with allophycocyanin (APC).

The three conjugates in TC690 have been produced from purified monoclonal mouse antibodies of isotypes IgG1, kappa (Anti-CD3), IgG1, kappa (Anti-CD19) and IgG1, kappa (Anti-CD45). TC690 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 50 tests (20 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).

Specificity

Anti-CD3, UCHT1, was included in the First and Third International Workshops and Conferences on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD3 (5). Anti-CD3, UCHT1, reacts with the ε-chain of CD3 (6).

Anti-CD19, HD37, was included in the Second, Third, Fourth and Fifth International Workshops and Conferences on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD19 (7, 8). Anti-CD19, HD37, labels human B cells in peripheral blood, bone marrow and other tissues. Additionally, B-cell lymphoproliferative disorders gave positive reactions with the antibody, i.e. acute lymphoblastic leukaemia, chronic lymphocytic leukaemia, hairy cell leukaemia, lymphoblastic lymphoma (Burkitt type), centroblastic/centrocytic lymphoma, and diffuse non-Hodgkin's lymphoma (7, 8). Anti-CD19, HD37, inhibits anti-Ig-induced B-cell activation and proliferation (2).

Anti-CD45, 2D1, was included in the Third International Workshop and Conference on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD45 (9).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.

Staining procedure

1. Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
2. Add 20 µL of TC690 and mix gently by using a vortex mixer.
3. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.
4. Add 100 µL of Uti-Lyse™ Reagent A (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 1 mL of Uti-Lyse™ Reagent B (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
7. Add 2 mL of PBS (Dako code No. S3024) to the tube and resuspend the cells by using a vortex mixer.
8. Repeat step 6.
9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.
10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the conjugated antibody reagent. The recommended Triple-Colour Control Reagent for TC690 is Dako code No. X0978.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, e.g. Dako EasyLyse™, code No. S2364, the PBS in step 9 should contain 1% paraformaldehyde or be analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis. The combination of FITC, RPE and APC is especially easy to work with since the emission spectra of these fluorochromes show minimum spectral overlap.

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1. Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematolymphoid neoplasias: results of an international consensus meeting. *Cytometry* 2001;46:23-7.
2. Pezzutto A, Dörken B, Rabinovitch PS, Ledbetter JA, Moldenhauer G, Clark EA. CD19 monoclonal antibody HD37 inhibits anti-immunoglobulin-induced B cell activation and proliferation. *J Immunol* 1987;138:2793-9.
3. Sato S, Tedder TF. BC3. CD19 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference*; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 133-5.
4. Leong AS-Y, Cooper K, Leong FJW-M. *Manual of diagnostic antibodies for immunohistology*. London: Greenwich Medical Media Ltd.; 2003. p. 121-4.
5. McMichael AJ, Gotch FM. T-cell antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. *Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference*; 1986 Sept 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 31-62.
6. Tunncliffe A, Olsson C, Traunecker A, Krissansen GW, Karjalainen K, de la Hera A. T3.2. The majority of CD3 epitopes are conferred by the epsilon chain. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. *Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference*; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 295-6.
7. Nadler LM. B cell/leukemia panel workshop: summary and comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. *Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leukocyte Differentiation Antigens*; 1984 Sept 17-20; Boston, USA. New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 3-43.
8. Mason DY, Ladyman H, Gatter KC. Immunohistochemical analysis of monoclonal anti-B cell antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. *Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leukocyte Differentiation Antigens*; 1984 Sept 17-20; Boston, USA. New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 245-55.
9. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. *Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference*; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 788-803.

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Agilent Technologies Singapore (International) Pte Ltd.