**MultiMix™**  
**Triple-Colour Control Reagent**  
**Mouse IgG1/FITC**  
**Rabbit F(ab')<sub>2</sub>/RPE**  
**Rabbit F(ab')<sub>2</sub>/APC**  
**Code No./ Réf./ Code-Nr. X0979**

**ENGLISH**

<b>Intended use</b>	For in vitro diagnostic use. X0979 is intended for use in flow cytometry. X0979 has been produced especially as a flow cytometry control reagent for Dako triple-colour reagents consisting of monoclonal mouse antibody of isotype IgG1 conjugated with FITC, and F(ab') <sub>2</sub> fragmented polyclonal rabbit antibodies conjugated with RPE and APC.
<b>Reagent provided</b>	X0979 comprises the following three matched fluorescent reagents: Monoclonal mouse IgG1 conjugated with fluorescein isothiocyanate isomer 1 (FITC). F(ab') <sub>2</sub> fragment of rabbit immunoglobulin conjugated with R-phycoerythrin (RPE). F(ab') <sub>2</sub> fragment of rabbit immunoglobulin conjugated with allophycocyanin (APC). Mouse IgG1/FITC is a purified monoclonal antibody, clone DAK-GO1, conjugated with FITC. The antibody specificity is directed towards <i>Aspergillus niger</i> glucose oxidase, an enzyme which is neither present nor inducible in mammalian tissues. Rabbit F(ab') <sub>2</sub> /RPE and Rabbit F(ab') <sub>2</sub> /APC have been prepared as follows: <ol style="list-style-type: none"> <li>1. The rabbit immunoglobulin (Ig) fraction is isolated from serum of healthy, non-immunized rabbits (Danish Whites) by salting out and ion exchange chromatography.</li> <li>2. Then the Ig fraction is solid-phase absorbed with human plasma. This reduces the non-specific background, and makes the reagent a better control for solid-phase absorbed antibodies.</li> <li>3. The immunoglobulin molecules are then degraded by pepsin and the F(ab')<sub>2</sub> fraction isolated by gel filtration.</li> <li>4. Finally, the F(ab')<sub>2</sub> fraction is conjugated with RPE or APC.</li> </ol> X0979 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2.
<b>Precautions</b>	<ol style="list-style-type: none"> <li>1. For professional users.</li> <li>2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>3. As with any product derived from biological sources, proper handling procedures should be used.</li> </ol>
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.
<b>Dilution guidelines</b>	In flow cytometry, X0979 may be used at a volume of 20 µL per test. This is a guideline only. The optimal amount may vary, e.g. depending on the specimen and preparation method, and should be determined by each individual laboratory. X0979 is an especially useful control reagent, since it aids in the discrimination between immunological and non-immunological reactivity of the test conjugate. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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Agilent Technologies Singapore (International) Pte Ltd.