Monoclonal Mouse Anti-Human CD19/APC, Clone HD37
Monoclonal Mouse Anti-Human CD19/FITC, Clone HD37
Monoclonal Mouse Anti-Human CD19/RPE, Clone HD37
Monoclonal Mouse Anti-Human CD19/RPE-Cy5, Clone HD37
Monoclonal Mouse Anti-Human CD19/PerCP-Cy5.5, Clone HD37

Code C7224
Code F0768
Code R0808
Code C7066
Code PR703

ENGLISH

Intended use

For in vitro diagnostic use.

C7224, F0768, R0808, C7066 and PR703 are intended for use in flow cytometry. CD19 is the broadest lineage-specific surface marker for B cells and it is present on the surface of virtually all B lymphocytes, including early B progenitor cells (1). CD19 expression is maintained in B-lineage cells that have undergone neoplastic transformation (2). Antibodies to CD19 are considered essential for the initial evaluation of acute and chronic lymphoproliferative disorders (3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Summary and explanation

CD19 is a 120 kDa transmembrane glycoprotein with a reduced Mr of 95 000 (2, 4). CD19 is a member of the immunoglobulin superfamily with two extracellular C2-type domains (2), and it is a critical signal transduction molecule that regulates B lymphocyte development, activation, and differentiation (4). CD19 expression is restricted to normal and neoplastic B cells, being absent from T cells, monocytes, and granulocytes (2). The CD19 antigen appears early during B cell maturation, probably at late pro-B cell stage around the time of Ig heavy chain rearrangement. It then persists during all stages of B cell maturation and is lost upon terminal differentiation to plasma cells (2). The advantage of immunodetection of CD19 expression is that B lineage leukemias and lymphomas rarely lose the epitope (5). Expression of CD19 in some acute myelocytic leukemia cells has been described (2).

Reagent provided

The Anti-CD19 conjugates, C7224, F0768, R0808, C7066 and PR703 have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains conjugate for 100 tests of up to 10⁶ leucocytes from normal human peripheral blood. See table 1 for volume of conjugate per test.

Isotype: IgG1, kappa. Conjugate concentration: See label on vial.

Table 1

Antibody Dako Code	Fluorochrome	Volume of conjugate per test	Isotype Reagent Dako Code
F0768	FITC (Fluorescein Isothiocyanate Isomer 1)	10 µL	X0927
R0808	RPE (R-Phycoerythrin)	10 µL	X0928
C7066	RPE-Cy5 (R-Phycoerythrin-Cyanine 5)	10 µL	X0955
C7224	APC (Allophycocyanin)	10 µL	X0968
PR703	PerCP-Cy5.5 (peridinin chlorophyll protein-Cyanine 5.5)	5 µL	No Dako Code*

*It is recommended to use a non-reactive monoclonal antibody of the same isotype and conjugated to PerCP-Cy5.5

Immunogen

Hairy cell leukemia cells (6).

Specificity

Anti-CD19, HD37, was included in the Second, Third, Fourth and Fifth International Workshops and Conferences on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD19 (7, 8).

Anti-CD19, HD37, labels human B cells in peripheral blood, bone marrow and other tissues. Additionally, B-cell lymphoproliferative disorders gave positive reactions with the antibody, i.e. acute lymphoblastic leukemia, chronic lymphocytic leukemia, hairy cell leukemia, lymphoblastic lymphoma (Burkitt type), centroblastic/centrocytic lymphoma, and diffuse non-Hodgkin's lymphoma. The antibody was shown to be unreactive with other cells in the human hematopoietic system and did not react with non-hematopoietic cells, e.g. in kidney, liver, breast or lung tissues (8). Anti-CD19, HD37, inhibits anti-Ig-induced B-cell activation and proliferation (1).

Precautions

1. For in vitro diagnostic use.
2. For professional users.
3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
4. As with any product derived from biological sources, proper handling procedures should be used.
5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
6. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Staining procedure

1. Collect venous blood into a test tube containing an anticoagulant.
2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
4. Mix 100 µL cell suspension with the volume of fluorochrome-conjugated Anti-CD19 stated in table 1.
5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
6. Incubate in the dark at 4 °C for 30 minutes.
7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Product-specific limitations

It has been observed that RPE-Cy5-conjugates and PerCP-Cy5.5-conjugates may bind to monocytes, resulting in background staining (9, 10).

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

1. Pezzutto A, Dörken B, Rabinovitch PS, Ledbetter JA, Moldenhauer G, Clark EA. CD19 monoclonal antibody HD37 inhibits anti-immunoglobulin-induced B cell activation and proliferation. *J Immunol* 1987;138:2793-9.
2. Sato S, Tedder TF. BC3. CD19 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan.* New York, London: Garland Publishing Inc.; 1997. p. 133-5.
3. Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematolymphoid neoplasias: results of an international consensus meeting. *Cytometry* 2001;46:23-7.
4. Sato S, Tedder TF. CD guide. CD19. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan.* New York, London: Garland Publishing Inc.; 1997. p. 764-5.
5. Leong AS-Y, Cooper K, Leong FJW-M. *Manual of diagnostic antibodies for immunohistology.* London: Oxford University Press; 1999. p. 67-8.
6. Pezzutto A, Dörken B, Feller A, Moldenhauer G, Schwartz R, Wernet P, et al. HD37 monoclonal antibody: a useful reagent for further characterization of "non-T, non-B" lymphoid malignancies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. *Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leukocyte Differentiation Antigens; 1984 Sept 17-20; Boston, USA.* New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 391-402.
7. Nadler LM. B cell/leukemia panel workshop: summary and comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. *Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leukocyte Differentiation Antigens; 1984 Sept 17-20; Boston, USA.* New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 3-43.
8. Mason DY, Ladyman H, Gatter KC. Immunohistochemical analysis of monoclonal anti-B cell antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, editors. *Leucocyte typing II. Proceedings of the 2nd International Workshop on Human Leukocyte Differentiation Antigens; 1984 Sept 17-20; Boston, USA.* New York, Berlin, Heidelberg, Tokyo: Springer-Verlag; 1986. Volume 2. p. 245-55.
9. van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). *Blood* 1996;88:2358-61.
10. Shapiro HM. *Practical Flow Cytometry,* John Wiley and Sons Inc.; 2003. 4th ed. p. 337.

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Agilent Technologies Singapore (International) Pte Ltd.