**Monoclonal Mouse
Anti-Human
CD22/PerCP-Cy5.5
Clone 4KB128**

Code PR707

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human CD22/PerCP-Cy5.5, Clone 4KB128, is intended for use in flow cytometry. Anti-CD22, clone 4KB128, is intended for use in the identification of cells expressing CD22. CD22 is a pan-B cell marker enabling detection of B cells in peripheral blood (1). CD22 is a recommended marker for B-cell lineage acute leukemia (2) together with a panel of other antibodies (3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.</p>
Synonyms for antigen	<p>B lymphocyte cell adhesion molecule (BL-CAM), Lyb-8.2 (1, 4).</p>
Summary and explanation	<p>CD22 is a 140 kDa single chain type 1 transmembrane glycoprotein (5). Cloning of the CD22 gene has revealed that CD22 is a member of the IgG superfamily and homologous to several other proteins, including myelin basic protein and members of the carcinoembryonic antigen (CEA) family (4). The structural features of CD22 suggest that CD22 would mediate intercellular adhesion, however, further work is needed to understand the in vivo function of CD22 and to identify physiologically significant ligands (1).</p> <p>CD22 expression is restricted to normal and neoplastic B cells and is absent from other hematopoietic cell types (2). CD22 is first expressed in pro-B and pre-B cells, largely as a cytoplasmic protein, and only later at the mature B-cell stage is it found on the cell surface (1, 2). The antigen is lost during the terminal stages of differentiation prior to the plasma cell stage. Most precursor-B-cell acute lymphoblastic leukemias (ALL) express both cytoplasmic and membrane forms of CD22, thus prominent surface expression of CD22 was observed in 90% of 52 precursor-B-ALL (2).</p>
Reagent provided	<p>The Anti-CD22 conjugate, PR707, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (5 µL conjugate /test) of up to 10⁶ leucocytes from normal human peripheral blood.</p> <p><u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.</p> <p>It is recommended to use a non-reactive monoclonal antibody of the same isotype and conjugated to PerCP-Cy5.5 as a negative control reagent.</p>
Immunogen	<p>Neoplastic cells from hairy cell leukemia (6).</p>
Specificity	<p>Anti-CD22, 4KB128, was included in the Third and Fifth International Workshops and Conferences on Human Leucocyte Differentiation Antigens (2, 6), and studies by different laboratories confirmed its reactivity with the CD22 antigen.</p> <p>Anti-CD22, 4KB128, was shown to detect surface membrane CD22 on precursor B cells (2). A weak staining of neutrophils in blood smears may occasionally be observed with this antibody.</p> <p>Similarly, a weak staining of neutrophils may be seen when clone 4KB128 is used in flow cytometry for intracellular staining of leucocytes from normal human peripheral blood.</p>
Precautions	<ol style="list-style-type: none">1. For in vitro diagnostic use.2. For professional users.3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.4. As with any product derived from biological sources, proper handling procedures should be used.5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	<p>Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.</p>
Staining procedure	<ol style="list-style-type: none">1. Transfer 100 µL of anticoagulated blood to a 12 mm x 75 mm polystyrene test tube.2. Add 5 µL of fluorochrome-conjugated anti-CD22 and mix gently by using a vortex mixer.3. Incubate the tube in the dark at 2-8 °C for 30 minutes.4. Add 2 mL of EasyLyse™ lysing reagent (Code S2364, diluted 20 times) to the tube and immediately mix gently using a vortex mixer. Incubate the tube for 15 minutes at room temperature in the dark.5. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.6. Add 2 mL of PBS (Code S3024) to the tube and resuspend the cells using a vortex mixer.

7. Repeat step 5.
8. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.4 mL PBS.
9. Analyze the sample on a flow cytometer.

Procedural notes

The recommended anticoagulant is EDTA. Heparin and citrate may be used but occasionally give more debris. It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Product-specific limitations

It has been observed that RPE-Cy5-conjugates and PerCP-Cy5.5-conjugates may bind to monocytes, resulting in background staining (7, 8).

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

References/ Bibliographie/ Literaturnachweise

- Tedder TF, Tuscano J, Sato S, Kehrl JH. CD22, a B lymphocyte-specific adhesion molecule that regulates antigen receptor signaling. *Annu Rev Immunol* 1997;15:481-504.
- Kehrl JH. B6. CD22 workshop panel report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference*; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press. Volume 1. 1995. p. 523-5.
- Wood, B. L., Arroz, M., Barnett, D., DiGiuseppe, J., Greig, B., Kussick, S. J., Oldaker, T., Shenkin, M., Stone, E. and Wallace, P. (2007), 2006 Bethesda International Consensus recommendations on the immunophenotypic analysis of hematolymphoid neoplasia by flow cytometry: Optimal reagents and reporting for the flow cytometric diagnosis of hematopoietic neoplasia. *Cytometry*, 72B: S14-S22.
- Wilson GL, Fox CH, Fauci AS, Kehrl JH. cDNA cloning of the B cell membrane protein CD22: a mediator of B-B cell Interactions. *J Exp Med* 1991;173:137-46.
- Engel P. CD guide. CD22. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. *Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference*; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 767-8. (3)
- Moldenhauer G, Schwartz R, Dörken B, Hämmerling GJ. B2.2. Biochemical characterization and epitope analysis of B-lymphocyte-specific surface antigens defined by clustering Workshop monoclonal antibodies. In: McMichael AJ, Beverley PCL, Cobbold S, Crompton MJ, Gilks W, Gotch FM, et al., editors. *Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference*; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 378-82.
- van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). *Blood* 1996;88:2358-61.
- Shapiro HM. *Practical Flow Cytometry*, John Wiley and Sons Inc.; 2003. 4th ed. p. 337.

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Agilent Technologies Singapore (International) Pte Ltd.