**Monoclonal Mouse
Anti-Human
CD1a/PerCP-Cy5.5
Clone NA1/34**

Code PR710

ENGLISH

Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human CD1a/PerCP-Cy5.5, Clone NA1/34, is intended for use in flow cytometry. Anti-CD1a, clone NA 1/34, is intended for use in the identification of cells expressing CD1a. Antibodies to CD1a may be a useful tool for identification of Langerhans' cell histiocytosis and for the classification of thymomas and malignancies of T-cell precursors (1) together with a panel of other antibodies (2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
Synonyms for antigen	R4 (3).
Summary and explanation	The human CD1 locus contains five genes, CD1A-E. CD1 molecules consist of glycosylated polypeptide α chains that span the cell membrane and are non-covalent associated with β 2-microglobulin (4). Anti-CD1a recognizes a heavy chain molecule of 49 kDa. This antigen is expressed strongly on cortical thymocytes and on a variety of antigen-presenting cells including Langerhans' cells and interdigitating dendritic cells. CD1a is not expressed on early thymocytes and is absent on mature peripheral blood T cells but intracytoplasmic expression is detected on activated T lymphocytes (1, 3, 4).
Reagent provided	The Anti-CD1a conjugate, PR710, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN_3 , pH 7.2. Each vial contains 100 tests (5 μL conjugate /test) of up to 10^6 leucocytes from normal human peripheral blood. <u>Isotype:</u> IgG2a, kappa. <u>Conjugate concentration mg/L:</u> See label on vial. It is recommended to use a non-reactive monoclonal antibody of the same isotype and conjugated to PerCP-Cy5.5 as a negative control reagent.
Immunogen	Human thymocytes (5).
Specificity	Anti-CD1a, NA1/34, was included in the First, Second, Third and Fourth International Workshop and Conferences on Human Leucocyte Differentiation Antigens, and studies by different laboratories confirmed its reactivity with the CD1a antigen (3).
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN_3), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Staining procedure	<ol style="list-style-type: none"> 1. Transfer 100 μL of anticoagulated blood to a 12 mm x 75 mm polystyrene test tube. 2. Add 5 μL of fluorochrome-conjugated anti-CD1a and mix gently by using a vortex mixer. 3. Incubate the tube in the dark at 2-8 °C for 30 minutes. 4. Add 2 mL of EasyLyse™ lysing reagent (Code S2364, diluted 20 times) to the tube and immediately mix gently using a vortex mixer. Incubate the tube for 15 minutes at room temperature in the dark. 5. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 μL of fluid in the tube. 6. Add 2 mL of PBS (Code S3024) to the tube and resuspend the cells using a vortex mixer. 7. Repeat step 5. 8. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.4 mL PBS. 9. Analyse the sample on a flow cytometer.

Procedural notes	The recommended anticoagulant is EDTA. Heparin and citrate may be used but occasionally give more debris. It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.
Product-specific limitations	It has been observed that RPE-Cy5-conjugates and PerCP-Cy5.5-conjugates may bind to monocytes, resulting in background staining (6, 7).

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
References/ Bibliographie/ Literaturnachweise

1. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 2003. p. 59-60.
2. Wood, B. L., Arroz, M., Barnett, D., DiGiuseppe, J., Greig, B., Kussick, S. J., Oldaker, T., Shenkin, M., Stone, E. and Wallace, P. (2007), 2006 Bethesda International Consensus recommendations on the immunophenotypic analysis of hematolymphoid neoplasia by flow cytometry: Optimal reagents and reporting for the flow cytometric diagnosis of hematopoietic neoplasia. *Cytometry*, 72B: S14–S22
3. Fainboim L, Salamone MC. CD guide. CD1a-e. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. *Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom.* New York: Oxford University Press Inc.; 2002. p. 747-8.
4. Fainboim L, Salamone MC. TC2. CD1 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan.* New York, London: Garland Publishing Inc.; 1997. p. 33-7.
5. McMichael AJ, Pilch JR, Galfré G, Mason DY, Fabre JW, Milstein C. A human thymocyte antigen defined by a hybrid myeloma monoclonal antibody. *Eur J Immunol* 1979;9:205-10.
6. van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). *Blood* 1996;88:2358-61.
7. Shapiro HM. *Practical Flow Cytometry*, John Wiley and Sons Inc.; 2003. 4th ed. p. 337.

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Agilent Technologies Singapore (International) Pte Ltd.