**Monoclonal Mouse
Anti-Human
CD7/PerCP-Cy5.5
Clone CBC.37**

Code PR711

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human CD7/PerCP-Cy5.5, Clone CBC.37, is intended for use in flow cytometry. Anti-CD7, Clone CBC.37, is intended for use in the identification of cells expressing CD7. Anti-CD7 is considered essential for the initial evaluation of T cell acute lymphoblastic leukemias (T-ALL) and T cell chronic leukemias (1) together with a panel of other antibodies (2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.</p>
Synonyms for antigen	gp40 (3).
Summary and explanation	<p>CD7 is a membrane-bound 40 kDa glycoprotein expressed on thymocytes, mature T cells, a large majority of natural killer cells, pluripotent hematopoietic stem cells, and progenitor cells of lymphoid and myeloid cells (3-6). CD7 is one of the earliest T-cell-specific antigens to be expressed in lymphocytes, and the only early marker to persist throughout differentiation (5).</p> <p>CD7 is expressed on most immature, common thymocytic and mature T-ALL, T-cell prolymphocytic leukemias, adult T-cell leukemias, large granular leukemias and some cutaneous T-cell leukemias/lymphomas (1, 7). CD7 is also expressed on special cases of acute myeloid leukemia of hematopoietic stem cell origin (8).</p> <p>CD7-deficient T cells have been observed in severe combined immunodeficiency (SCID) and rheumatoid arthritis (3, 7).</p>
Reagent provided	<p>The Anti-CD7 conjugate, PR711, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (5 µL conjugate /test) of up to 10⁶ leucocytes from normal human peripheral blood.</p> <p><u>Isotype:</u> IgG2b, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.</p> <p>It is recommended to use a non-reactive monoclonal antibody of the same isotype and conjugated to PerCP-Cy5.5 as a negative control reagent.</p>
Immunogen	CEM T-cell line (ATCC CCL-199), a T-lymphoblastoid cell line, established from a patient with acute lymphoblastic leukemia (9).
Specificity	<p>SDS-PAGE analysis of immunoprecipitates formed between CEM cell lysate and the antibody shows that the antibody precipitates a 40 kDa polypeptide corresponding to CD7 (9).</p> <p>In flow cytometry, the antibody strongly labels the cell surface of Jurkat and CEM T-cell lines, whereas the Raji B-cell line is negative. Additionally, the antibody blocks the reaction of the CD7 antibody, clone DK24 (9).</p>
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Staining procedure	<ol style="list-style-type: none"> 1. Transfer 100 µL of anticoagulated blood to a 12 mm x 75 mm polystyrene test tube. 2. Add 5 µL of fluorochrome-conjugated anti-CD7 and mix gently by using a vortex mixer. 3. Incubate the tube in the dark at 2-8 °C for 30 minutes. 4. Add 2 mL of EasyLyse™ lysing reagent (Code S2364, diluted 20 times) to the tube and immediately mix gently using a vortex mixer. Incubate the tube for 15 minutes at room temperature in the dark. 5. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube. 6. Add 2 mL of PBS (Code S3024) to the tube and resuspend the cells using a vortex mixer. 7. Repeat step 5. 8. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.4 mL PBS. 9. Analyse the sample on a flow cytometer.

Procedural notes The recommended anticoagulant is EDTA. Heparin and citrate may be used but occasionally give more debris. It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Product-specific limitations It has been observed that RPE-Cy5-conjugates and PerCP-Cy5.5-conjugates may bind to monocytes, resulting in background staining (10, 11).

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- Shapiro HM. Practical Flow Cytometry, John Wiley and Sons Inc.; 2003. 4th ed. p. 337.
Agilent Technologies Singapore (International) Pte Ltd.