**Polyclonal Rabbit  
Anti-Human  
Lambda Light Chains/PerCP-Cy5.5  
Rabbit F(ab')<sub>2</sub>**

**Code PR712**

**ENGLISH**

<b>Intended use</b>	For in vitro diagnostic use. PR712 is intended for use in flow cytometry. The antibody is intended for use in the identification of cells in peripheral blood that express lambda light chains. In flow cytometry, antibodies to lambda light chains are useful for the demonstration of cell surface lambda light chains, and thus, for the identification of monoclonality (clonal excess) in B-cell lymphoproliferative disorders when used together with a panel of other antibodies. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
<b>Summary and explanation</b>	Most B cells, except for pre-B progenitors, pre-B cells, and mature plasma cells, express immunoglobulin on their surface. Each cell expresses only one light chain type. In normal peripheral blood and lymph nodes, there is a mixture of kappa-positive and lambda-positive cells, with two-thirds of the cells expressing kappa and one-third expressing lambda (1). Since lymphoid neoplasms are usually clonal expansions of a single cell, malignant cells uniformly express the same light chain isotype. Neoplastic B-cell chronic lymphoproliferative disorders can frequently be suspected on basis of demonstration of a marked predominance of cells expressing a single light chain type (2, 3). Rare cases of kappa-negative or lambda-negative B-cell non-Hodgkin's lymphoma have been described (3, 4).
<b>Reagent provided</b>	The Anti-Lambda Light Chains conjugate, PR712, has been produced from a F(ab') <sub>2</sub> fragment of affinity-isolated polyclonal rabbit antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Use 5 µL of the product for staining of up to 10 <sup>6</sup> leucocytes from normal human peripheral blood. It is recommended to use a non-reactive polyclonal F(ab') <sub>2</sub> antibody conjugated to PerCP-Cy5.5 as a negative control reagent. <u>Conjugate concentration:</u> See label on vial.
<b>Manufacturing of the PerCP-Cy5.5 conjugate</b>	<ol style="list-style-type: none"> <li>1. The antibody used for PerCP-Cy5.5 conjugation has been solid-phase absorbed with human plasma proteins to remove traces of contaminating antibodies.</li> <li>2. The absorbed antibody has been further purified by affinity chromatography on a column with immobilized human lambda light chains.</li> <li>3. The affinity-isolated antibody has subsequently been degraded with pepsin and the F(ab')<sub>2</sub> fragment isolated by gel filtration.</li> <li>4. Finally, the F(ab')<sub>2</sub> fragment has been conjugated with tandem fluorochrome peridinin chlorophyll protein-Cyanine 5.5 (PerCP-Cy5.5).</li> </ol>
<b>Immunogen</b>	Polyclonal immunoglobulin chains of lambda type isolated from a pool of human sera.
<b>Specificity</b>	Anti-Lambda Light Chains reacts with free lambda chains as well as lambda chains in intact immunoglobulin molecules. The specificity of this product has been ascertained by crossed immune electrophoresis. To obtain maximum sensitivity, the crossed immune electrophoresis specificity test was performed prior to affinity purification and pepsin degradation. When Anti-Lambda Light Chains is applied as described in the staining procedure in combination with Anti-CD19/RPE, HD37 or Anti-CD19/FITC, HD37, on lysed human whole blood, a specific staining of a part of the CD19-positive B lymphocytes is seen corresponding to the expected range of lambda light chain expression. Flow cytometric analysis of single suspensions from formalin-fixed, paraffin-embedded tissue samples demonstrated that Anti-Lambda Light Chains labels reactive hyperplastic lymph nodes (10/10 cases) (1). In non-Hodgkin's lymphomas, the antibody labeled 5/10 cases. The remaining cases were positive for kappa light chains (4/10 cases) or showed no expression of light chains (1/10 case) (3). Flow cytometric analysis of peripheral blood lymphocytes demonstrated that Anti-Lambda Light Chains labels a proportion of B-cell chronic lymphocytic leukemias. Thus, in one study of 121 cases (5), 17 were positive for lambda. In another study of 165 cases (6), 68 were positive for lambda.
<b>Precautions</b>	<ol style="list-style-type: none"> <li>1. For in vitro diagnostic use.</li> <li>2. For professional users.</li> <li>3. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>4. As with any product derived from biological sources, proper handling procedures should be used.</li> <li>5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.</li> <li>6. Unused solution should be disposed of according to local, State and Federal regulations.</li> </ol>
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. During storage a small precipitate may occasionally develop causing a fine granular non-specific staining. By a simple filtration (0.22 µm cellulose acetate filter), the original high quality of the conjugate will be restored. Conjugates should not be stored in diluted form. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.
<b>Procedural notes</b>	Before staining samples of peripheral blood, the mononuclear cells must be isolated by centrifugation on a separation medium or the blood sample washed to remove the soluble serum proteins. As human monocytes bind serum immunoglobulins via their surface Fc receptors, these cells should be removed by depletion or identified. The recommended anticoagulant is EDTA. Heparin and citrate may be used but occasionally give more debris. It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis. It has been observed that PerCP-Cy5.5-conjugates may bind to monocytes, resulting in background staining (7, 8).
<b>Staining procedure</b>	<ol style="list-style-type: none"> <li>1. Collect venous blood into a test tube containing an anticoagulant.</li> <li>2. Transfer 100 µL of the anticoagulated blood into the test tube.</li> <li>3. Add 2 mL of 0.01 mol/L PBS, pH 7.4. Mix gently by using a vortex mixer.</li> </ol>

4. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 µL of fluid.
5. Repeat steps 3 and 4 two more times.
6. Add 10 µL rabbit immunoglobulin fraction, Code X0903, for blocking. Mix gently by using a vortex mixer and incubate in the dark at 37 °C for 30 minutes.
7. Add 5 µL of PerCP-Cy5.5-conjugated Anti-Lambda Light Chains. Mix gently.
8. Use a non-reactive F(ab')<sub>2</sub> fragment of rabbit immunoglobulin, and conjugated with the same fluorochrome, as a negative control.
9. Incubate in the dark at 4 °C for 30 minutes.
10. Add 1-2 mL of erythrocyte-lysing reagent to each tube and mix gently. Follow the reagent manufacturer's recommendations for volume, time and temperature of incubation. (During development product performance was validated with EasyLyse erythrocyte-lysing reagent Code S2364).
11. Centrifuge at 300 x g for 5 minutes.
12. Aspirate the supernatant, leaving approximately 50 µL of fluid.
13. Add 2 mL PBS, pH 7.4. Mix gently by using a vortex mixer.
14. Repeat steps 11 and 12.
15. Resuspend pellet in an appropriate fluid for flow cytometry, e.g. PBS, pH 7.4.
16. Analyze on a flow cytometer.

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
## References/ Bibliographie/ Literaturnachweise

1. Deegan MJ. B lymphocytes and plasma cells: their development and identification. In: Keren DF, editor. Flow cytometry in clinical diagnosis. Chicago: ASCP Press; 1989. p. 139-63.
2. Johnson A, Olofsson T. Flow cytometric clonal excess analysis of peripheral blood, routine handling, and pitfalls in interpretation. Cytometry 1993;14:188-95.
3. Leers MPG, Theunissen PHMH, Ramaekers FCS, Schutte B, Nap M. Clonality assessment of lymphoproliferative disorders by multiparameter flow cytometry of paraffin-embedded tissue: an additional diagnostic tool in surgical pathology. Hum Pathol 2000;31:422-7.
4. Kaleem Z, Zehnauer BA, White G, Zutter MM. Lack of expression of surface immunoglobulin light chains in B-cell non-Hodgkin lymphomas. Am J Clin Pathol 2000;113:399-405.
5. Cartron G, Linassier C, Bremond JL, Desablens B, Georget MT, Fimbel B, et al. CD5 negative B-cell chronic lymphocytic leukemia: clinical and biological features of 42 cases. Leuk Lymphoma 1998;31:209-16.
6. Gandini D, Lanza F, Latorraca A, Levato F, Del Senno L, Castoldi G. Immunophenotypic and genotypic characterization of B-cell chronic lymphocytic leukemia patients from Northern Italy. Haematologica 1993;78:18-24.
7. van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). Blood 1996;88:2358-61.
8. Shapiro HM. Practical Flow Cytometry, John Wiley and Sons Inc.; 2003. 4th ed. p. 337.

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Agilent Technologies Singapore (International) Pte Ltd.