**Monoclonal Mouse
Anti-Human
Plasma Cell/PerCP-Cy5.5**
Clone VS38c

Code PR713

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>PR713 is intended for use in flow cytometry for detection of plasma cells in peripheral blood. PR713 is intended to be used in a panel with other antibodies. The antibody clone VS38c reacts with normal and neoplastic human plasma cells. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.</p>
Synonym for antigen	<p>p63 (1, 2).</p>
Summary and explanation	<p>Anti-Human Plasma Cell recognizes the intracytoplasmic antigen p63 of 64 kDa present in normal and neoplastic cells (4). p63 is an intracellular type II transmembrane protein which is localized in the rough endoplasmic reticulum. The function of p63 is unknown, but its abundance in secretory cells and localization to the endoplasmic reticulum suggests a conserved role in protein processing or secretion (1).</p>
Reagent provided	<p>The Anti-Human Plasma Cell conjugate, PR713, has been produced from a purified monoclonal mouse antibody conjugated to peridinin chlorophyll protein-Cyanine 5.5 (PerCP-Cy5.5).</p> <p>The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (5 µL of conjugate per test for up to 10⁶ cells).</p> <p><u>Isotype</u>: IgG1, kappa.</p> <p><u>Conjugate concentration</u>: See label on vial.</p>
Immunogen	<p>MCF-7 AdR, a multi-drug resistant breast carcinoma cell line (3).</p>
Specificity	<p>Anti-Human Plasma Cell, VS38c, was submitted to and characterized at the sixth International Workshop and Conference on Human Leucocyte Differentiation Antigen (1). Anti-Human Plasma Cell, VS38c, labels normal and neoplastic plasma cells in various tissues, but does not label peripheral blood lymphocytes, although some weak labeling of monocytes and neutrophils could occasionally be detected (3). The VS38c antibody also labels melanocytic cells, particularly malignant melanoma cells, a number of epithelial cells from various organs (5) and osteoblasts from bone marrow biopsies (4). Myeloma, plasmacytoma, and lymphoplasmacytoid lymphomas are strongly labeled and about a third of large B-cell lymphomas are weakly labeled. T-cell lymphomas are negative. Occasional weak labeling of endothelial cells may occur (1).</p>
Precautions	<ol style="list-style-type: none">1. For in vitro diagnostic use.2. For professional users.3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.4. As with any product derived from biological sources, proper handling procedures should be used.5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	<p>Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.</p>
Procedural notes	<p>It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control.</p> <p>Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p> <p>It has been observed that PerCP-Cy5.5-conjugates may bind to monocytes, resulting in background staining (6, 7).</p>
Staining procedure	<ol style="list-style-type: none">1. Transfer 50 µL (up to 10⁶ cells) of the cell suspension to be analyzed (anticoagulated whole blood or mononuclear cells) to a test tube.2. Add 100 µL Dako IntraStain Reagent A (Fixation), Code K2311. Mix gently by using a vortex mixer to ensure that the cells are in suspension.

3. Incubate at room temperature for 15 minutes.
4. Add 2 mL PBS and mix gently by using a vortex mixer.
5. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 µL of fluid.
6. Mix thoroughly by using a vortex mixer to ensure that the cells are in suspension and add 100 µL Dako IntraStain Reagent B (Permeabilization), Code K2311. Add 5 µL of PR713. Mix gently by using a vortex mixer to ensure that the cells are in suspension.
7. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control.
8. Incubate in the dark at room temperature for 15 minutes.
9. Repeat steps 4 and 5.
10. Resuspend the pellet in an appropriate fluid for flow cytometric analysis, e.g. 0.3 mL PBS, pH 7.4.
11. Analyze on a flow cytometer

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7. Shapiro HM. Practical Flow Cytometry, John Wiley and Sons Inc.; 2003. 4th ed. p. 337.
Agilent Technologies Singapore (International) Pte Ltd.