



YAMASA CORPORATION

DIAGNOSTICS DEPARTMENT
Tokyo Office
1-23-8, Nihonbashi-kakigaracho
Chuo-ku, Tokyo, 103-0014, Japan
Tel. +81-3-3668-8558 Fax. +81-3-3668-8407
www.yamasa.com/shindan/index.htm

Rat Bone EDTA Extract

Code No. : 7609
Quantity : 0.5 ml (5 ng)
Quality : Extract from rat bone with EDTA solution.
Solution : 10 ng /ml in 15 mM sodium phosphate, 5 mM EDTA Na_4 ,
0.01% tween 20, 0.1% proclin, 150 mM NaCl, pH 6.8
Storage: The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.
Note : Research use only.